

THE ORIGINATOR®

Clear Aligner System



Care
Guide

Welcome to straight teeth without braces.



You have a lot to look forward to, congratulations on getting started with your new Originator Clear Aligner System. If you follow these simple care instructions, and wear your Originator as prescribed by your dental professional, you'll soon have the straight, beautiful smile you've always wanted.



How the Originator works.

The Originator is specifically designed to address minor corrections that don't require braces. It straightens your teeth by applying constant, gentle pressure to gradually move your teeth in stages, one Originator at a time. Each aligner repositions your teeth a bit farther, until you reach your goal. It's that simple. You'll be checking in with your dental professional on a regular basis to monitor your progress and make any necessary adjustments.



Wear your Originator as prescribed.

It's best to wear your Originator at all times, removing it only to eat, brush and floss. You may experience a slight discomfort at first, but this is normal and will quickly go away. You'll wear the same aligner as prescribed by your orthodontist before moving to the next tray in your series. Your dental professional has prescribed your wearing schedule, and you'll want to stick to it for the fastest possible results.



How to take care of your Originator.

Clean your Originator with a soft toothbrush, denture cleaner or mouthwash.

Toothpaste tends to scratch the surface and make it look “cloudy,” so try to avoid its use. Do not soak or place your aligner in hot water as the heat may melt or misshapen it. When not in your mouth, keep your Originator in its protective case and store at room temperature. Take extra caution never to leave it in your car or in the sun. Be sure to brush, floss and rinse regularly, ideally after every meal, to avoid discoloration, tooth decay or bad breath.



Trust the quality of your Originator.

You can be certain you’re wearing the very best clear aligners available today. They’re made of a durable medical-grade plastic featuring a smooth, anti-bacterial surface that keeps them extra clean. And they’re manufactured by TP Orthodontics, a leader in the orthodontic industry. Over the years, TPO has made more custom tooth positioning appliances than any other company.



TP Orthodontics, Inc.

Originator is a registered trademark of TP Orthodontics, Inc.

©2013 TP Orthodontics, Inc.

996-662ENG/SPN, Rev. 0