ENGLISH

Click-It® Forceps

The forceps frame is manufactured from stainless steel. The tips are made from advanced polymers. The forceps utilize replaceable tips to extend the life of the forceps while allowing optimum performance. Dry heat sterilization is recommended as autoclaving may result in rust, pits and fissures.

Tip Replacement

The tips should be inspected for wear with each use. They will need to be replaced as necessary to ensure proper function and to avoid damage to the Click-It brackets.

To remove a worn or flattened tip:

- Hold the forceps firmly in one hand while using a wedged blunt instrument such as a small flat blade screwdriver to pull up on the retaining tab until it breaks. **The broken tab should be contained and eye protection should be provided and worn.**
- After the tab is broken from the tip simply slide the worn tip off the forceps and discard.
- Hold the new tip between your thumb and forefinger and align the blade of the forceps into the opening. Slide the tip onto the blade until it is fully seated and the retaining tab snaps into the holes provided on the blade. Test the forceps for correct assembly. The tips on the forceps should align when the forceps are fully closed tip-to-tip. If they do not align, check each tip to be certain that it is secure and undamaged. Replace if necessary to avoid damage to the Click-It brackets.

Proper Care

For optimum care of this precision instrument, TPO recommends the following:

For Cleaning: Non-Ionic No-Rinse General Purpose Concentrate
For Sterilization: Rapid Dry Heat Sterilization System

For proper cleaning, removal of the tips from the forceps is not recommended.

Manual Cleaning

- Use a non-ionic no-rinse general purpose concentrate.
- Hand scrub forceps and tips to remove all fluids or other foreign substances from the forceps and the tips.
- Dry the forceps and tips thoroughly and remove any remaining liquid from hinge area.

Automated Cleaning

- Ultrasonically clean with no-rinse ultrasonic cleaning solution containing rust inhibitors.
- Clean the forceps and the tips for 5 minutes if using small ultrasonic cleaners and 10 minutes for large ultrasonic cleaners.
- Dry the forceps and the tips thoroughly and remove any remaining liquid from hinge area.

Sterilization

Autoclave

If autoclaving the forceps and tips add anti-rust tablets to water source and milk bath to lubricate.

If using a porous load autoclave or any other autoclave that creates a pre-vacuum the forceps and tips must also be properly wrapped before insertion into the autoclave.

Autoclave for 18 minutes at 121°C (249.8°F) or 3 minutes at 134°C-137°C (273.2°F-278.6°F) while avoiding any metal to metal contact.

After autoclaving, dry thoroughly and remove any remaining liquid from the hinge area and the tips (a blast of compressed air is an effective means of doing this).

Heat

TPO strongly recommends the use of Instrument Shield (Order Number 100-879) on the metallic forceps prior to heat sterilization.

TPO recommends dry heat.

For dry heat sterilization it is recommended that forceps and tips are heated to 150°C (302°F) for 150 minutes or 219°C (426°F) for 23 minutes for a rapid cycle dry heat sterilizer, brown spots may occur where foreign material is not completely removed before sterilization.

For best results, clean instruments before sterilizing.

DO NOT USE Iodophors (iodine) for cleaning prior to sterilization.

Immersion Disinfection

If using an immersion process to disinfect and/or sterilize instruments, please be aware of the following:

Specific products within a general group of chemical compounds (such as 2% glutaraldehydes) may be formulated quite differently, and will therefore provide very different results in terms of rust and corrosion. Use only products containing rust inhibitors and those which have been demonstrated as safe for use on carbon steel.

Also, keep in mind that chemical solutions which are diluted with water may experience alteration of their chemical balance by the additives in your local water supply.

- For best results, DO NOT immerse forceps and tips in these materials:
 - Glutaraldehydes with added phenols
 - Chlorine dioxide compounds
 - Alkali solutions
 - Acidic solution (below pH 6.0)
 - Highly reactive fluoridated compound
- After rinsing disinfectant from instruments, dry thoroughly and remove any remaining liquid from the hinge areas and the tips (a blast of compressed air is an effective means of doing this).
- After rinsing disinfectant from instruments, thoroughly dry tips and remove moisture from hinge areas (a blast of compressed air is an effective means of doing this).

There is no limit to reprocessing and end of life is determined by wear and damage caused by use.

To properly maintain forceps and tips do the following:

- Clean and sterilize forceps and tips as soon as it is practical after use and store forceps in a clean dry area.
- Periodically examine your instruments for any signs of unusual wear or discoloration:
 - Check for damage or wear on replaceable tips
 - Check for any signs of corrosion

NOTICE:

TP Orthodontics, Inc. warrants that its products are free from defects in materials and workmanship. TPO's liability and the customer's sole remedy in the event of any claimed defect shall be limited to replacement of the product, or at TPO's discretion, refund of the purchase price. This warranty is in lieu of any other warranties, express or implied, including any implied warranty of merchantability, fitness for a particular purpose, or otherwise, and of any other obligations or liability on TPO's part. Other than the warranties listed above, under no circumstances will TPO be liable for any loss, damages arising in connection with use, or inability to use, TPO's products.

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