

INSTRUCTIONS

Series 151

tportho.com

CrossLink is a registered trademark of TP Orthodontics, Inc.

© 2006 TP Orthodontics, Inc.
996-274, Rev. 1

ENGLISH

Instructions

1. Thoroughly clean teeth to be banded with pumice using a rotating bristle brush or rubber cup. Rinse teeth thoroughly with water, isolate teeth and air dry.
2. Acid conditioning of the enamel is not required for cementing bands.
3. Optional: For maximum bond strength use CrossLink Cement Intensifier (151-282). Dispense a few drops on a mixing pad and apply Intensifier to the inside of the band with a small bristle brush.
4. Apply a small bead of cement to the inside of the band. Spread the paste with a small spatula to evenly coat the entire band interior.

NOTE: CrossLink Cement is light sensitive. If several bands are filled at the same time, cover bands to protect from direct light.
5. Fully seat band and clean off all excess flash.
6. Light cure by placing light tip close to occlusal surface of tooth and moving light slowly around at the top of the band. Cure for 30 seconds. To assure proper curing, the light should be able to provide a minimum of 400 milliwatts of continuous output.

Helpful Tips:

CrossLink will also bond to plastic and composite restorations when those surfaces have been conditioned with Light Cure Adhesion Promoter.

Storage:

Best temperature for use is 68°-77°F (20°-25°C). Storage temperature is best at 35°-80°F (2°-27°C). CrossLink has a 3-year shelf life without refrigeration.

Note: This is a single-use device. Any reuse of this product will run the risk of cross-contamination to patients.

NOTICE:

TP Orthodontics, Inc. warrants that its products are free from defects in materials and workmanship. TPO's liability and the customer's sole remedy in the event of any claimed defect shall be limited to replacement of the product, or at TPO's discretion, refund of the purchase price. This warranty is in lieu of any other warranties, express or implied, including any implied warranty of merchantability, fitness for a particular purpose, or otherwise, and of any other obligations or liability on TPO's part. Other than the warranties listed above, under no circumstances will TPO be liable for any loss, damages arising in connection with use, or inability to use, TPO's products.

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