INSTRUCTIONS

Preformed Palatal Arch Bar

Series 352

ENGLISH

The removable Preformed Palatal Arch Bar is designed to rotate, expand or torque maxillary first molars. These movements can be accomplished individually, dually, or all three simultaneously.

To use this appliance, it is necessary to attach a molar sheath on the lingual of each maxillary first molar band. Molar sheaths measuring 0.91mm x 2.03mm (.036" x .080") I.D. should be attached at the same occluso-gingival and mesio-distal position as the buccal tubes.

Measure

1. With a flexible millimeter ruler, measure from the middle of one sheath to the middle of the opposite sheath. Select the appropriate Palatal Arch Bar.

Form

2. If not enough palatal clearance, select the smaller size bar and expand slightly. If too much clearance, select the next larger bar and slightly constrict.

Activate for Expansion

3. Palatal Arch Bars are activated for expansion by opening the loop slightly. Although the loop is compressed while seating the bar in the molar sheaths, it will return to its activated position and expand the molars bilaterally.

NOTE: Expanding the loop will also cause the double-back ends of the bar to toe-out at the distal, causing the molars to rotate. If rotation is not desired, adjust the distal ends of the double-back mesially.

Activate for Rotation

4. To activate Palatal Arch Bars for rotation, simply adjust the distal end of the double-back lingually or buccally before engagement in the molar sheaths.

Activate for Torque

5. The double-back ends of the Palatal Arch Bar make it possible to torque the molar roots lingually toward the palate or buccally toward the buccal plate. Simply grasp the double-back ends with ribbon arch pliers or universal pliers and incorporate lingual or buccal torque.

Placement


6. Seat one side of the Palatal Arch Bar into the lingual molar sheath and check for desired activation before seating the other side. Activation may include rotation and expansion for increasing arch length or torque for increasing molar anchorage. These forces can be applied individually or simultaneously.

7. Engage the remaining side to the molar sheath.

NOTE: The Palatal Arch Bar should be securely ligated into the lingual sheath with steel or elastic ligatures.

8. Adjust the Palatal Arch Bar with the same frequency as a maxillary archwire or facebow if either is present.

Mechanical Force Illustration

NOTE: This is a single-use device. Any reuse of this product will run the risk of cross-contamination to patients.

Warning: This stainless steel product contains nickel and chromium. Direct patient to consult physician if allergic reaction occurs.

NOTICE:

TP Orthodontics, Inc. warrants that its products are free from defects in materials and workmanship. TPO's liability and the customer's sole remedy in the event of any claimed defect shall be limited to replacement of the product, or at TPO's discretion, refund of the purchase price. This warranty is in lieu of any other warranties, express or implied, including any implied warranty of merchantability, fitness for a particular purpose, or otherwise, and of any other obligations or liability on TPO's part. Other than the warranties listed above, under no circumstances will TPO be liable for any loss, damages arising in connection with use, or inability to use, TPO's products.

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