INSTRUCTIONS

M-D Reducer Mesiodistal Tooth Reduction Instrument

Series 100

Symbols Used on Labeling									
Manufacturer:	Date of manufacture:	Authorized Representative in the European Community:	Batch code:	Caution! See warnings and precautions:	Consult instructions for use:	Do not re-use:	Use by date:	Do not use if package is damaged:	Product is not sterilized:

TP Orthodontics, Inc.
100 Center Plaza
La Porte, Indiana 46350-9672 USA

tportho.com

Advena Ltd.
Tower Business Centre, 2nd Flr.
Tower Street, Swatar, BKR 4013
Malta

M-D Reducer developed by Dr. Harry Barrer.

© 2006 TP Orthodontics, Inc.

996-532, Rev. C - 4/21

ENGLISH

The M-D Reducer is a clinically proven method to remove interproximal tooth mass in flat or contoured planes. It simulates the normal attritional process by reducing teeth mesiodistally to produce additional intra-arch space. The M-D Reducer may also be used to open contact points, recontour restorations or remove overhanging fillings. The M-D Reducer consists of a frame with interchangeable aluminum oxide cutter blades. The cutter blades fasten on two simple frame pins for controlled flat or curved cutting strokes.


Instructions for Use

- For tight contacts, hand-hold frame and pull blade through contact point to create enough space to position cutter without jamming.
- To cut flat planes, hold handle steady in horizontal plane. Do not flex blade against interproximal surfaces.
- To recontour, press laterally on handle to flex blade around interproximal of tooth and gently rock horizontally to produce a rounded surface.
- The same procedures can be used for anterior or posterior teeth.
- Site of enamel reduction can be controlled by using single- or double-sided blades.
- Coarse blades are normally used for cutting while fine blades are used for polishing or for opening up tight contacts.
- Avoid striking opposing incisors with frame. Work with patient's mouth wide open, and protect moving end of frame from hard or soft tissue with finger.
- Make full inciso-gingival cuts—do not create gingival shelves.

Mesiodistal Reduction Techniques

Reduce tooth mass and create "right-angle" interproximal surfaces.

"Keystoning"—create angulated, flat interproximal surfaces to resist tendency for relapse after orthodontic treatment.

The M-D Reducer Offers Many Advantages:

- Simple and rapid action.
- Complete control of cutting blade.
- Complete shape and contour control.
- Safe—minimizes soft tissue damage.
- Excellent patient acceptance.
- Choice of blades—abrasive one side or both sides—coarse for cutting and fine for polishing.
- Blades will cut amalgam, silicates or composites.
- Use of two blades at same time doubles thickness of cut.

Sterilization Instructions

To sterilize your M-D Reducer handle, use cold sterilization.

NOTE: M-D Reducer blades are recommended for single use only. Any reuse of this product will run the risk of cross-contamination to patients.

NOTICE:

TP Orthodontics, Inc. warrants that its products are free from defects in materials and workmanship. TPO's liability and the customer's sole remedy in the event of any claimed defect shall be limited to replacement of the product, or at TPO's discretion, refund of the purchase price. This warranty is in lieu of any other warranties, express or implied, including any implied warranty of merchantability, fitness for a particular purpose, or otherwise, and of any other obligations or liability on TPO's part. Other than the warranties listed above, under no circumstances will TPO be liable for any loss, damages arising in connection with use, or inability to use, TPO's products.

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