
INSTRUCTIONS

Reverse Palatal Arch Bar

Series 352

Symbols Used on Labeling

| | |
|--|--|
| | Manufacturer. |
| | Date of manufacture. |
| | Batch code. |
| | Caution! See warnings and precautions. |
| | Consult instructions for use. |
| | Do not re-use. |
| | Do not use if package is damaged. |
| | Product is not sterilized. |

TP Orthodontics, Inc.
100 Center Plaza
La Porte, Indiana 46350-9672 USA

tportho.com

ENGLISH

The removable Palatal Arch Bar is designed to rotate, expand or torque maxillary first molars. These measurements can be accomplished individually, dually, or all three simultaneously.

To use this appliance, it is necessary to attach a molar sheath on the lingual of each maxillary first molar band. Molar sheaths measuring .036" x .72" I.D. should be attached at the same occluso-gingival and mesio-distal position as the buccal tubes.

Measure

- With a flexible millimeter ruler, measure from the middle of one sheath to the middle of the opposite sheath. Follow the contour of the palate without impinging on the palatal tissue. Allow 1.5 to 2 mm clearance between ruler and palate. Select the appropriate Palatal Arch Bar.

Form

- Grasp the bar at the vertical loop with pliers that are wide enough to cover the entire loop. Push down on one side of the bar, then the other to form a slight "V" shape.
- Contour the Palatal Arch Bar to the shape of the palate with contouring pliers. Leave 1.5 to 2 mm clearance between bar and palate.

Activate for Expansion

- Palatal Arch Bars are activated for expansion by opening the loop slightly. Although the loop is compressed while seating the bar in the molar sheaths, it will return to its activated position and expand the molars bilaterally.

NOTE: Expanding the loop will also cause the double-back ends of the bar to toe-in at the distal, causing the molars to rotate. If rotation is not desired, adjust the distal ends of the double-back buccally.

Activate for Rotation

- To activate Palatal Arch Bars for rotation, simply adjust the distal end of the double-back lingually or buccally before engagement in the molar sheaths.

Activate for Torque

- The double-back ends of the Palatal Arch Bar make it possible to torque the molar roots lingually toward the palate or buccally toward the buccal plate. Simply grasp the double-back ends with ribbon arch pliers or universal pliers and incorporate lingual or buccal torque.

Placement


- Seat one side of the Palatal Arch Bar into the lingual molar sheath and check for desired activation before seating the other side. Activation may include rotation and expansion for increasing arch length or torque for increasing molar anchorage. These forces can be applied individually or simultaneously.

8. Engage the remaining side into the molar sheath.

WARNING: The Palatal Arch Bar should be securely ligated into the lingual sheath with steel or elastic ligatures.

9. Adjust the Palatal Arch Bar with the same frequency as a maxillary archwire or facebow if either is present.

Mechanical Force Illustration

NOTE: This is a single-use device. Any reuse of this product will run the risk of cross-contamination to patients.

NOTICE:

TP Orthodontics, Inc. warrants that its products are free from defects in materials and workmanship. TPO's liability and the customer's sole remedy in the event of any claimed defect shall be limited to replacement of the product, or at TPO's discretion, refund of the purchase price. This warranty is in lieu of any other warranties, express or implied, including any implied warranty of merchantability, fitness for a particular purpose, or otherwise, and of any other obligations or liability on TPO's part. Other than the warranties listed above, under no circumstances will TPO be liable for any loss, damages arising in connection with use, or inability to use, TPO's products.

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