# INSTRUCTIONS

## Ceramic Brackets

Click the language you wish to view for Ceramic Brackets to be taken directly to that page.

- English
- Dutch
- French
- German
- Italian
- Portuguese
- Spanish
- Chinese
- Japanese

TP Orthodontics, Inc.  
100 Center Plaza  
La Porte, Indiana 46350-9672  
USA

[tportho.com](http://tportho.com)


0297 Advena Ltd.  
Tower Business Centre, 2nd Flr.  
Tower Street, Swatar, BKR 4013 Malta

## Ceramic Brackets

The products covered by this instruction include:

- ClearVu® Ceramic Brackets
- InVu® Ceramic Brackets
- Tip-Edge® PLUS Ceramic Brackets

### Symbols Used on Labeling

	Manufacturer.
	Date of manufacture.
	Authorized Representative in the European Community.
	Batch code.
	Caution! See warnings and precautions.
	Consult instructions for use.
	Do not re-use.
	Use by date.
	Do not use if package is damaged.
	Medical device.
	Prescription device.

 **TP Orthodontics, Inc.**  
100 Center Plaza  
La Porte, Indiana 46350-9672  
USA

[tportho.com](http://tportho.com)

 0297 Advena Ltd.

Tower Business Centre, 2nd Flr.  
Tower Street, Swatar, BKR 4013 Malta

ClearVu, InVu, and Tip-Edge are registered trademarks of TP Orthodontics, Inc. Patents apply.

### Indications for Use

The Ceramic Bracket is a component of the fixed orthodontic appliance, which is one of the treatment forms available for correcting or improving malocclusions (mal positioning of the teeth). The ceramic orthodontic bracket functions as a point for attachment for the orthodontic arch wire, lock pins, or any necessary auxiliaries onto the tooth. The bracket is bonded to the tooth and they serve primarily as a mean of transferring applied stresses to achieve a clinically desired movement of the teeth.

### Warnings

- This product is intended for single use. Its reuse may cause cross contamination, loss of material integrity, and compromise overall product safety and performance.
- Don't reuse the bracket in the event of a premature debond as the bracket may have been deformed. Bond a new bracket (with the same part number).
- Ceramic is such a hard material that abrasion of the enamel can occur if incisal surfaces occlude. Each patient's occlusion should be taken into consideration when determining if and when ceramic can be safely bonded and still avoid contact with opposing teeth.
- Do not touch polymer mesh bases prior to bonding.
- Remove any stray adhesive from bracket surface and markings prior to curing. Cured adhesive covering any bracket identification or positional marking will prevent removal of these markings.
- Do not force auxiliaries or archwires into ceramic brackets.
- Steel ligatures should not be tied as tightly on ceramic as on stainless steel.
- Do not use adhesion promoters or enhancers on enamel.
- The color mark provides orientation reference for the positioning of the bracket.

### Precautions

Advise patients not to bite hard food or materials, to prevent damage to the brackets and/or teeth.

### Contraindications

- Bonding brackets to patients with compromised oral health, such as poor oral hygiene, periodontal disease, or preexisting dental conditions (i.e. decalcification, carious lesions, large restorations) can increase the risk of tooth damage.
- Bonding brackets to uncooperative patients can increase the risk of tooth damage.

### Instructions for Use

This product must be used by qualified orthodontists or dentists trained in orthodontics.

### BRACKET BONDING

**Please note:** Do not use adhesion promoters or enhancers on enamel.

No special bonding technique or bonding system is required for polymer mesh bases. Follow the adhesive manufacturer's instructions for a traditional mechanical bond, just like metal brackets. The patented polymer mesh base has undercuts that enable a mechanical interlocking to occur.

Polymer mesh bases can be safely bonded and debonded from crowns and other non-enamel surfaces. No special conditioners or primers are required. When bonding to gold, amalgam, stainless steel, acrylic or composite restorations, use of an adhesion promoter is acceptable.

The tooth surface should be cleaned before bonding the bracket. Adhesive manufacturer's instructions should be followed.

Remove any adhesive flash to minimize adhesive staining, plaque retention, and to facilitate debonding.

After the bracket is bonded, remove the positioning jig (for brackets that have a jig).

## Ceramic Brackets

The products covered by this instruction include:

- ClearVu® Ceramic Brackets
- InVu® Ceramic Brackets
- Tip-Edge® PLUS Ceramic Brackets

### *Positional Markings*

Ceramic brackets may include temporary markings within the horizontal and vertical slots that serve as positional guides during bonding. These markings may be removed after bonding, and prior to installing the archwire, by cleaning with a standard toothbrush, toothpaste and water.

### *Rebonding*


In the case of a premature debond, do not rebond. The polymer mesh base would become damaged when debonding and would not bond to the enamel properly if reused.

Clean excess adhesive from tooth, prepare the tooth surface according to the adhesive manufacturer's instructions, and bond a new bracket.


### **BRACKET DEBONDING**

Polymer mesh bases debond without the need for special instruments for techniques. The polymer mesh safety base allows safe and easy debonding. The debonding instrument of choice is a ligature cutting plier. Squeeze the mesiodistal edges of the mesh base, grasping as closely as possible to the enamel surface. Do not grasp the bracket itself. The polymer mesh base will easily release from the tooth while retaining most of the adhesive.

As an added safety feature, it is very helpful to leave brackets ligated to the archwire when debonding. If excessive bonding adhesive prevents applying pressure to the base or tends to deflect the beaks of the ligature cutting pliers away from the polymer mesh base, remove the adhesive from the mesiodistal surfaces of the base with a steel, fissure bur. Avoid using a carbide bur or diamond instrument as they are not necessary and could accidentally mar the enamel surface. Once the adhesive is removed, apply compressive pressure directly to the edges of the base near the enamel.

Apply slow, steady compressive force to the mesiodistal edges of the bracket's mesh base, making sure to be close to the enamel. The base will deform and debond.

Care should be taken to contain small parts so they are not accidentally inhaled or ingested by the patient. It is helpful to leave these ligated to the archwire when debonding.

### **Storage Condition**

Product should be kept in a clean location.

### **Shelf Life**

5 years.

### **REPORTING SERIOUS INCIDENTS**

User and Patient Notice: If a serious incident occurs with this medical device, report it to the manufacturer and to the competent medical authority for the user/patient's country.

#### **NOTICE:**

TP Orthodontics, Inc. warrants that its products are free from defects in materials and workmanship. TPO's liability and the customer's sole remedy in the event of any claimed defect shall be limited to replacement of the product, or at TPO's discretion, refund of the purchase price. This warranty is in lieu of any other warranties, express or implied, including any implied warranty of merchantability, fitness for a particular purpose, or otherwise, and of any other obligations or liability on TPO's part. Other than the warranties listed above, under no circumstances will TPO be liable for any loss, damages arising in connection with use, or inability to use, TPO's products.

## Ceramic Brackets

The products covered by this instruction include:

- ClearVu® Ceramic Brackets
- InVu® Ceramic Brackets
- Tip-Edge® PLUS Ceramic Brackets

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