

The Flip-Lock[®] HERBST[®]

Only from TP Orthodontics

The Flip-Lock HERBST from TP Orthodontics is quickly and easily activated using the patented "Flip-Lock" mechanism. There are no pins or screws. Activation is simple — just place crimpable spacers at regular intervals onto the exposed piston.

The Flip-Lock HERBST is available with and without a cantilever arm. Optional archwire slots allow the appliance to be placed in conjunction with brackets.




Component Parts PACKAGE OF 10

- 610-403** Pistons
- 610-404** Tubes (5 left, 5 right)
- 610-405** Ball Hinge with small pad WITHOUT SLOT
- 610-422** Posterior Ball Hinge WITHOUT SLOT and large welding pad
- 610-423** Posterior Ball Hinge WITH SLOT and large welding pad
- 610-421** Cantilever/anterior ball assembly WITHOUT SLOT (5 left, 5 right)
- 610-425** One-piece Cantilever/anterior ball assembly WITH SLOT (5 left 5 right)



Crimpable Spacers PACKAGE OF 10

- 610-401** 1 mm
- 610-402** 2 mm
- 610-406** 3 mm
- 610-407** 4 mm
- 610-408** 5 mm



Order online at
www.tportho.com

Single Patient and Inventory Kits

- 610-400** Single Patient Appliance Kit (2 pistons, 2 tubes, 4 ball hinges on small pad WITHOUT SLOT)
- 610-410** Single Patient Appliance Kit (2 pistons, 2 tubes, 4 ball hinges on large pad WITH SLOT)
- 610-430** Single Patient Appliance Kit (2 pistons, 2 tubes, 4 ball hinges on large pad WITHOUT SLOT)
- 610-420** Single Patient Appliance Kit (2 tubes, 2 pistons, 2 posterior ball hinges, 2 cantilever/ball hinges WITHOUT SLOT)
- 610-424** Single Patient Appliance Kit (2 tubes, 2 pistons, 2 posterior ball hinges, 2 cantilever/ball hinges WITH SLOT)
- 610-409** Inventory Kit (30 pistons, 30 tubes, 60 ball hinges, 95 spacers)

