



## NV<sup>®</sup> Hybrid Tubes

**are the only low-profile molar tubes with a polymer mesh base.**

The NV Hybrid tube features a translucent polymer mesh base that allows light to penetrate completely when curing.

This results in a stronger bond which increases stability and reliability throughout the course of treatment.

# Single

## 6/6 Upper First Molar

Archwire Tube Size	Tube Length	Torque	Tip	Distal Offset	Left/Right	Bondable	
						Pack of 10	Pack of 4
0.46 x 0.64mm (.018" x .025")	4.0	-10°	0°	6°	L	<b>391-912</b>	<b>391-4912</b>
					R	<b>391-913</b>	<b>391-4913</b>
	-14°	0°	10°	L	<b>391-916</b>	<b>391-4916</b>	
				R	<b>391-917</b>	<b>391-4917</b>	
0.56 x 0.71mm (.022" x .028")	4.0	-10°	0°	6°	L	<b>391-914</b>	<b>391-4914</b>
					R	<b>391-915</b>	<b>391-4915</b>
	-14°	0°	10°	L	<b>391-918</b>	<b>391-4918</b>	
				R	<b>391-919</b>	<b>391-4919</b>	

Package of 10

## 6/6 Lower First Molar

Archwire Tube Size	Tube Length	Torque	Tip	Distal Offset	Left/Right	Bondable	
						Pack of 10	Pack of 4
0.46 x 0.64mm (.018" x .025")	4.0	-20°	0°	0°	L	<b>391-920</b>	<b>391-4920</b>
					R	<b>391-921</b>	<b>391-4921</b>
	-25°	0°	0°	L	<b>391-924</b>	<b>391-4924</b>	
				R	<b>391-925</b>	<b>391-4925</b>	
0.56 x 0.71mm (.022" x .028")	4.0	-20°	0°	0°	L	<b>391-922</b>	<b>391-4922</b>
					R	<b>391-923</b>	<b>391-4923</b>
	-25°	0°	0°	L	<b>391-926</b>	<b>391-4926</b>	
				R	<b>391-927</b>	<b>391-4927</b>	

Package of 10

**To order, visit [tportho.com](http://tportho.com) or call your Customer Service Representative.**

**Australia:** 1 800 643 055

**United Kingdom:** 44 (0) 113 2539192

**South Africa:** +27 (0)11 100 0956

**United States:** 800-348-8856

NV is a registered trademark of TP Orthodontics, Inc. and manufactured under US Patents 5,263,859, 6,685,468, 6,746,242, 7,780,442, 8,925,719; Britain 1545380, 1562511; China 03822189.6, ZL 03824889.1; Europe 1545380, 1562511; France 1545380, 1562511; Germany 1545380, 60332263.8-08; Italy 1545380, 1562511; Japan 4414334, 4503440; Mexico 265725, 274623; Spain ES2362089T3, 1562511.

© 2010 TP Orthodontics, Inc.

POD\_NV Hybrid Tubes\_ENG\_2015, Rev. 0



**TP Orthodontics, Inc.**