

Nu-Edge® Self-Ligating Mini Cobalt Chromium Brackets

MBT* McLaughlin, Bennett, Trevisi System

PATIENT KITS

Slot Size	With Hooks	
0.46mm (.018")	Cuspid/Bicuspid	293-291SL
0.56mm (.022")	Cuspid/Bicuspid	293-294SL

Kits contain 20 brackets.



MAXILLARY BRACKETS

	Slot Size	Left/Right	Bondable
Central	0.46mm (.018")	L	293-310SL
Torque +17°		R	293-311SL
Angulation + 4°	0.56mm (.022")	L	293-312SL
		R	293-313SL
Lateral	0.46mm (.018")	L	293-314SL
Torque +10°		R	293-315SL
Angulation + 8°	0.56mm (.022")	L	293-316SL
		R	293-317SL
Cuspid with Distal Hook	0.46mm (.018")	L	293-322SL
Torque -7°		R	293-323SL
Angulation +8°	0.56mm (.022")	L	293-324SL
		R	293-325SL
Bicuspid with Distal Hook	0.46mm (.018")	L	293-330SL
Torque - 7°		R	293-331SL
Angulation 0°	0.56mm (.022")	L	293-332SL
		R	293-333SL

Package of 5


MANDIBULAR BRACKETS

	Slot Size	Left/Right	Bondable
Anteriors	0.46mm (.018")	L/R	293-354SL
Torque -6°			
Angulation 0°	0.56mm (.022")	L/R	293-355SL
Cuspid with Distal Hook	0.46mm (.018")	L	293-360SL
Torque -6°		R	293-361SL
Angulation +3°	0.56mm (.022")	L	293-362SL
		R	293-363SL
1st Bicuspid with Distal Hook	0.46mm (.018")	L	293-268SL
Torque -12°		R	293-269SL
Angulation +2°	0.56mm (.022")	L	293-270SL
		R	293-271SL
2nd Bicuspid with Distal Hook	0.46mm (.018")	L	293-276SL
Torque -17°		R	293-277SL
Angulation +2°	0.56mm (.022")	L	293-278SL
		R	293-279SL

Package of 5

nu-edge® SL

Cobalt Chromium Brackets



Bracket Identification Chart

*MBT is a registered trademark of 3M Unitek. TP Orthodontics version of the McLaughlin, Bennett and Trevisi prescription is not intended to be a duplication of, nor an endorsement by, these doctors and such endorsement is not implied by the use of this term.

Nu-Edge® Self-Ligating Cobalt Chromium Brackets

Compatible with the Roth system: not endorsed by Dr. Roth

PATIENT KITS

Slot Size	With Hooks	
0.46mm (.018")	Cuspid/Bicuspid	293-106SL
0.56mm (.022")	Cuspid/Bicuspid	293-107SL

Kits contain 20 brackets.



MAXILLARY BRACKETS

		Slot Size	Left/ Right	Bondable
Central		0.46mm	L	293-010SL
	Torque	+17°	R	293-011SL
	Tip	+5°	L	293-012SL
		0.56mm	R	293-013SL
Lateral		0.46mm	L	293-014SL
	Torque	+8°	R	293-015SL
	Tip	+9°	L	293-016SL
		0.56mm	R	293-017SL
Cuspid with Distal Hook		0.46mm	L	293-022SL
	Torque	0°	R	293-023SL
	Tip	+11°	L	293-024SL
		0.56mm	R	293-025SL
Bicuspid with Distal Hook		0.46mm	L	293-098SL
	Torque	-7°	R	293-099SL
	Tip	0°	L	293-100SL
		0.56mm	R	293-101SL

Package of 5

MANDIBULAR BRACKETS

		Slot Size	Left/ Right	Bondable
Anteriors		0.46mm	L/R	293-050SL
	Torque	0°		
	Tip	0°	L/R	293-052SL
		0.56mm		
Cuspid with Distal Hook		0.46mm	L	293-058SL
	Torque	-11°	R	293-059SL
	Tip	+5°	L	293-060SL
		0.56mm	R	293-061SL
1st Bicuspid with Distal Hook		0.46mm	L	293-090SL
	Torque	-17°	R	293-091SL
	Tip	0°	L	293-092SL
		0.56mm	R	293-093SL
2nd Bicuspid with Distal Hook		0.46mm	L	293-082SL
	Torque	-22°	R	293-083SL
	Tip	0°	L	293-084SL
		0.56mm	R	293-085SL

Package of 5