

Submitted

Digital Impression Model
 Upper Lower

Remove all attachments from model:


Yes No

TRAYS

Upper Lower
 .030** .040"
 *.030" is the standard thickness.


STRIP TEETH

None
 Strip 0. ___ mm from each tooth at contact points marked.



RESET TEETH

None Over rotate
 Compromise Ideal
 Reset teeth circled Overcorrect labiolingually



VIEW AND APPROVE DIGITAL SET UP

Yes – view final only
 Yes – all steps
 No



SPECIAL INSTRUCTIONS

Originator and TPO are registered trademarks of TP Orthodontics, Inc.

ORIGINATOR® Rx
Aligner Tray System

The Originator is used to correct minor to moderate anterior crowding or spacing. It can also be used to correct minor orthodontic relapse.

The Originator Basic System includes up to 10 trays—including 5 dual arches (upper and lower). Most minor to moderate corrections can be achieved in three trays.



Right ←————→ Left

View your Originator set-up online at
www.tportho.com/custom or
call TPO® for assistance at 800-348-8856.

Packaging Instructions

PVS Impressions

1. Wrap each impression in protective packaging material such as bubble wrap or foam.
2. Place each half in special divided shipping carton with completed prescription form and ship to TPO.

Models

1. Wrap each half of the model in the special polyurethane square TPO provides and secure with a rubber band.
2. Place each half in the special divided shipping carton with prescription form and ship to TPO.

Caution: Place only one set of models in each shipping carton. Be sure that divider separates each half.

Mailing & Shipping
 100 Center Plaza
 La Porte, Indiana 46350-9672 USA
 Phone: 800-348-8856
 219-785-2591

Instructions for digital/online
prescription submissions
can be found at

tportho.com/custom

This is my first case with TP Orthodontics.

ACCOUNT NO. _____

Dr. _____

Address _____

City _____

State _____

Zip _____ Phone () _____

E-Mail (required) _____

Patient's Name _____

Date Shipped to TPO _____

Date Required _____

This custom-made device is manufactured to satisfy the design characteristics and properties specified by the prescribing doctor for this specific patient, and is intended for the exclusive use of the named patient.

PLEASE DO NOT WRITE IN THIS SPACE			