

# Pre-Finisher<sup>®</sup>

## Preformed Finishing Appliance

### Efficient finishing for aligner cases.

The Pre-Finisher appliance is a semi-custom solution perfect for improving the occlusion in aligner and fixed-appliance patients. When used as an adjunct to aligner treatment, they can reduce or eliminate additional trays required for case refinement – while correcting the posterior open bite often seen with aligner patients.



Nothing could be easier! Near the end of treatment, take a simple measurement and select a Pre-Finisher appliance from your in-office inventory kit – no impressions required. You'll appreciate the time-savings and the final results.

- No laboratory wait or fees; immediate delivery to the patient
- Relates teeth in both arches simultaneously; corrects open bite
- Aligns, rotates and closes spaces
- Sized in one-millimeter increments for precise fitting
- Optional seating springs facilitate placement and retention
- Ideal solution for the already-compliant aligner patient
- Results in 2-3 weeks

### Real results:

#### Case A

Before Pre-Finisher delivery



After 3 weeks of night time wear



(Improved Class II and open bite)

#### Case B

Before Pre-Finisher delivery



(Improved open bite)

After 3 weeks of wear



The Refine<sup>®</sup> Orthodontic Finishing Appliance is also available as a customized solution specially developed for aligner treatment.

Find out more at [tportho.com](http://tportho.com)  
or call **800-348-8856**



# Pre-Finisher® Preformed Finishing Appliance

## Parts List

### Multi-Appliance Kits

#### Convincer Kits

(10 Assorted Sizes)	Plain	With Seating Springs
Nonextraction	190-010	190-110
Four Bicuspid Extraction	190-410	190-510
Two Upper Bicuspid Extraction	190-210	190-310

#### Standard Kits

(15 Assorted Sizes)	Plain	With Seating Springs
Nonextraction	190-015	190-115
Four Bicuspid Extraction	190-415	190-515
Two Upper Bicuspid Extraction	190-215	190-315

#### Thrift Kits

(20 Assorted Sizes)	Plain	With Seating Springs
Nonextraction	190-020	190-120
Four Bicuspid Extraction	190-420	190-520
Two Upper Bicuspid Extraction	190-220	190-320



Size, in millimeters (last two digits of part number) molded into each appliance.

### Individual Appliances

The last two digits of the order number indicate size of Pre-Finisher appliance in millimeters.

#### Nonextraction Series

Can also be used on first molar extraction cases.

Plain					With Seating Springs				
190-040	190-044	190-048	190-052	190-056	190-140	190-144	190-148	190-152	190-156
190-041	190-045	190-049	190-053	190-057	190-141	190-145	190-149	190-153	190-157
190-042	190-046	190-050	190-054	190-058	190-142	190-146	190-150	190-154	190-158
190-043	190-047	190-051	190-055		190-143	190-147	190-151	190-155	

Nineteen sizes (40 mm to 58 mm).

#### Four Bicuspid Extraction Series

Plain					With Seating Springs				
190-444	190-448	190-452	190-456	190-460	190-544	190-548	190-552	190-556	190-560
190-445	190-449	190-453	190-457		190-545	190-549	190-553	190-557	
190-446	190-450	190-454	190-458		190-546	190-550	190-554	190-558	
190-447	190-451	190-455	190-459		190-547	190-551	190-555	190-559	

Seventeen sizes (44 mm to 60 mm).

#### Two Upper Bicuspid Extraction Series

Plain					With Seating Springs				
190-244	190-247	190-250	190-253	190-256	190-344	190-347	190-350	190-353	190-356
190-245	190-248	190-251	190-254		190-345	190-348	190-351	190-354	
190-246	190-249	190-252	190-255		190-346	190-349	190-352	190-355	

Thirteen sizes (44 mm to 56 mm).

### Asian Appliances

Asian appliances are specially designed to address common Asian dental anatomy traits:

- Increased incisal mamelon ridge definition
- Decreased clinical crown length

#### Asian Nonextraction Series

Plain			
190-848	190-850	190-852	190-854
190-849	190-851	190-853	190-855

Eight sizes (48 mm to 55 mm).

#### Asian Four Bicuspid Extraction Series

Plain				
190-644	190-647	190-650	190-653	190-656
190-645	190-648	190-651	190-654	190-657
190-646	190-649	190-652	190-655	

Fourteen sizes (44 mm to 57 mm).