

eXact[®]

Clear Orthodontic Adhesive

eXact Clear Orthodontic Adhesive is a medium-high viscosity light cure adhesive. This innovative adhesive system can be used for metal, plastic or ceramic attachments and cures completely clear - making it an exceptional system for aesthetic brackets.

eXact Clear Orthodontic Adhesive has exceptional workability/tack properties to prevent bracket drift, produces little flash and is extremely easy to clean up.



eXact Clear Light Cure Kit

151-270

Kit Contains:

- Light Cure Syringes, 6 g (.212 oz.), 3
- Replacement Tips
- TPO[®] Blue Etchant Gel in Syringe, 5 cc
- Light Cure Sealant, 6 cc
- Disposable Brushes, 50
- Foam Pellets, 100
- Spatula, 1
- Work Pad, 50 sheets

Replacement Materials

- 151-271** Paste in Syringe, 6 g (.212 oz.)
- 150-713** TPO Blue Etchant Gel in Syringe, 5 cc
- 151-272** Light Cure Sealant, 6 cc
- 150-704** Disposable Brushes, 50
- 150-705** Foam Pellets, 150
- 150-904** White Spatulas, 50
- 150-903** Work Pad, 50 sheets
- 100-620** Adhesive Mixing and Orientation Trays, 25

Features

- Excellent tack properties to prevent bracket sliding
- Cures clear and produces little flash for minimal clean up
- Complete wetting of the bracket mesh and prepared tooth surface for a secure bond
- For use on metal or ceramic brackets

Find out more at tportho.com
or call **800-348-8856**




CrossLink[®]

Light Cure Band Cement

CrossLink light cure band cement is available in no-mix, no-mess, single-use dispensing caplets. Disposable squeeze caplet has an angled, tapered stem for excellent control when dispensing paste. Simply cut the caplet tip and dispense cement into bands. For faster dispensing, the tip can be cut shorter to produce a larger diameter opening.

Fluoride ion releasing paste does not require sealant, and there is no rush to beat setting time. Complete cure is accomplished in 30 seconds. CrossLink provides excellent band retention.

CrossLink is available in a natural or blue color. Blue helps clinicians to see cement when cleaning excess flash after cementing bands.



Natural Color Cement

- 151-287** Large Patient Caplets, 12 caplets, 1.8 g (.063 oz.) each (4 molar bands)
- 151-283** Small Patient Caplets, 30 caplets, .4 g (.014 oz.) each (1 molar band)
- 151-280** CrossLink Kit - Natural, 24 g (.847 oz.) paste

Blue Color Cement

- 151-288** Large Patient Caplets, 12 caplets, 1.8 g (.063 oz.) each (4 molar bands)
- 151-284** Small Patient Caplets, 30 caplets, .4 g (.014 oz.) each (1 molar band)
- 151-285** CrossLink Kit - Blue, 24 g (.847 oz.) paste

Replacement Materials

- 151-281** Natural Paste in Syringe, 4 g (.141 oz.)
- 151-286** Blue Paste in Syringe, 4 g (.141 oz.)
- 150-904** White Spatulas, 50

Single Patient Caplets

- No-mix, no-mess, single-use dispensing caplets
- Disposable squeeze caplet with angled, tapered stem for excellent control

Preloaded Syringes

- No sealant required
- Complete cure in 30 seconds
- Excellent band retention
- Available in natural or blue



Find out more at tportho.com
or call **800-348-8856**



TP Orthodontics, Inc.
tportho.com