


The future of Class II correction is here.

The Jasper Vektor is the **ONLY** appliance that utilizes gentle, intrusive force vectors – just what is needed for efficient Class II correction.



- Produces gentle, curved force vectors since the appliance flexes in 25-45% of its overall length

- Robust nickel titanium design

- Easy installation

- Free from the “food zone”

- Allows natural lateral mandibular movement for patient comfort

Other appliances produce extrusive tipping forces on the maxilla and anterior teeth, which is always contraindicated for the over erupted Class II patient. Since the vector control module (VCM) is made from nickel titanium, the appliance is break-resistant and force values of only 3.5 ounces are produced. Provide your patients the latest Class II correction alternative with the Jasper Vektor appliance.


See why the Jasper Vektor appliance is in a class of its own. Call or visit us online:

800-348-8856 | www.tportho.com/vektor

Jasper Vektor® Appliance

The Jasper Vektor Class II Correction Appliance is a fixed intraoral device that achieves class II correction easily and efficiently with low forces. The patented design of the appliance introduces gentle intrusive vector forces - unlike other class II correctors that place extrusive vector forces on the anterior teeth.

The appliance produces gentle sweeping forces, and the mandible is free to move laterally resulting in patient comfort. The appliance is easy to install, clears the “food zone” and is compliance-independent.



Patented

Kits	Order No.
Class II Corrector Multi-Kit Kits consists of 1 pair size 27, 2 pairs size 30, 2 pairs size 33, 1 pair size 36, 1 archwire installation tool, 20 archwire locks	610-500
Class II Corrector Kit—size 27	610-527
Class II Corrector Kit—size 30	610-530
Class II Corrector Kit—size 33	610-533
Class II Corrector Kit—size 36	610-536

Single Kit includes Left / Right Component plus 4 Archwire Locks



Multi-Kit

Component Part	Order No.
27mm Left	610-527L
27mm Right	610-527R
30mm Left	610-530L
30mm Right	610-530R
33mm Left	610-533L
33mm Right	610-533R
36mm Left	610-536L
36mm Right	610-536R
Package of 5	
Archwire Locks	610-502
Package of 10	

Tool	Order No.
Archwire Lock Installation Tool	610-501
Each	



Left / Right Component



Archwire Lock



Archwire Lock
Installation Tool



TP Orthodontics, Inc.
tportho.com