


# CROSSLINK<sup>®</sup> BAND CEMENT PATIENT CAPLETS

ONLY FROM  
TP ORTHODONTICS



TP Orthodontics, Inc.

# CROSSLINK<sup>®</sup> BAND CEMENT

- No Mix
- No Measure
- No Waste
- No Cross-Contamination
- No Rush
- No Wait

**Time Saving and Convenient**



# CROSSLINK® BAND CEMENT

TPO is the first to offer light cure band cement packaged in single-patient caplets.

## **Time Saving**

CrossLink band cement helps reduce chair time. This light-activated, glass ionomer hybrid cement requires no mixing and no measuring. CrossLink band cement cures in 30 quick seconds, and requires no conditioners or sealants.

CrossLink also contains fluoride to control decalcification under the bands during treatment.

## **Convenient**

CrossLink band cement is pure convenience. The caplets eliminate excessive waste and cross contamination between patients.

After each patient, the caplet is thrown away.

With single-patient CrossLink caplets, simply cut the tip and squeeze. The narrow-tipped design accurately dispenses adhesive directly into the molar bands.

The flow of cement can be customized by cutting the tip shorter to produce a larger-diameter opening.



### Two Caplet Sizes

CrossLink caplets are designed for dispensing cement into a single patient's band or bands. Small CrossLink caplets are designed to fill a single molar band, while the large size contains cement sufficient for four bands.



Large CrossLink caplet -  
1.8g each



Small CrossLink caplet -  
.4g each

### Two Color Choices

CrossLink band cement is available in two colors. Natural cement blends with enamel, while blue cement helps expedite clean-up after band placement.



Natural-colored cement



Blue-colored cement

At TP Orthodontics we never stop looking for a better way to do things. Orthodontists around the world choose TPO brands because they provide innovative orthodontic solutions that no other manufacturer can match.

Our engineers understand the difference between technology that works in a lab and what really works in your practice. This means not a single product is released into the market until we know it will perform to your highest expectations.

We believe in developing products that make your treatment more effective and your patients happier. Our dedicated customer service personnel are standing by to assist you. When you choose TPO as your orthodontic partner, you get only the best for your practice.

## ONLY THE BEST FOR YOUR PRACTICE

ONLY FROM TP ORTHODONTICS

**Order these and all TPO products  
online at [www.tportho.com](http://www.tportho.com) or call**

USA 800-348-8856

England 44 (0) 113 2539192

Australia 1-800 643 055 or 61 +3 9342 3200

South Africa +27 (0) 11 510 0097



TP Orthodontics, Inc.

*Headquarters*

100 Center Plaza

La Porte, Indiana 46350-9672 USA

ONLY FROM  
**TP ORTHODONTICS**