

SmileTracker™
MOBILE PATIENT MONITORING

- A dynamic approach to orthodontic treatment planning
- More timely adjustment of appliances when required
- Interactive and engaging for patients
- Secure and compliant

Set your practice apart with **SmileTracker**

Smile-Tracker.com



An Interactive Process for Creating Better Smiles

SmileTracker is a convenient way for you to monitor your patients' progress and compliance in between appointments helping you keep your treatment plans on target. **SmileTracker's** interactive platform allows orthodontists to view current patient photos and up-to-date time-lapse video footage of their treatment progress. Patients will stay engaged throughout treatment with a fun app that rewards them, and even lets them share their orthodontic treatment via social media. Use of an Apple iPhone is required.




Doctor logs into *Smile-Tracker.com*, creates a new account and issues an access code to the patient.



Patient downloads iPhone app, enters access code and begins tracking progress by taking daily photos.



Doctor logs into *Smile-Tracker.com*, views patient photos and a time-lapse video of the patient's progress.



SmileTracker™

MOBILE PATIENT MONITORING

Sign up for a **60-day FREE trial** and keep track of your patients between appointments with **SmileTracker!**



Orthodontist Pricing*

5 Users

\$30
per month

10 Users

\$60
per month

25 Users

\$150
per month

50 Users

\$300
per month

75 Users

\$450
per month

100 Users

\$600
per month

*Access to your patients is free for 60 days, then auto-billed to your credit card. See website for details.

tportho.com/smiletracker



Brought to you by

TP Orthodontics, Inc.