

1

AESTHETIC BRACKET SYSTEMS



TP Orthodontics, Inc.

tportho.com

cosmeticbraces.com

AESTHETIC BRACKET SYSTEMS

ClearVu®

- Mini, Roth
- Mini, MBT

ClearVu® with Pre-Applied Adhesive

- Mini, Roth
- Mini, MBT

InVu®

- Mini, Roth
- Roth
- Mini, MBT
- MBT
- Standard Edgewise

InVu® with Pre-Applied Adhesive

- Mini, Roth
- Roth
- Mini, MBT
- MBT

Click-It® Self-Ligating

- MBT

Tip-Edge PLUS®

Tip-Edge PLUS® with Pre-Applied Adhesive



Personalized
Color-Matching
Technology®

TPO Aesthetic Brackets

feature *Personalized Color-Matching Technology®* which allows our brackets to blend with your patient's tooth color.

ClearVu® Mini Ceramic Brackets

Compatible with the Roth system: not endorsed by Dr. Roth

PATIENT KITS

	With Hooks	0.46mm (.018") Slot	0.56mm (.022") Slot
Kits contain all ceramic brackets			
Upper 3-3	Cuspid	297-100CV	297-101CV
Upper 5-5	Cuspid/Bicuspid	297-404CV	297-405CV
Upper & Lower 3-3	Cuspid	297-102CV	297-103CV
Upper & Lower 5-5	Cuspid/Bicuspid	297-406CV	297-407CV
Kits combine ceramic with metal* brackets			
Up. 5-5 Ceramic with Lo. 5-5 Metal	Cuspid/Bicuspid	297-051CV	297-056CV

*Nu-Edge Cobalt Chromium brackets.

297-490CV Bi-Dimensional Brackets with Hooks
 .018" centrals and laterals with .022" cuspids and bicuspids
 Kit contains 20 brackets.



Patented

MAXILLARY BRACKETS

	Slot Size	Left/Right	Bondable
Central	0.46mm (.018")	L	297-138CV
Torque +14°		R	297-139CV
Tip +5°	0.56mm (.022")	L	297-150CV
		R	297-151CV
Lateral	0.46mm (.018")	L	297-148CV
Torque +7°		R	297-149CV
Tip +8°	0.56mm (.022")	L	297-152CV
		R	297-153CV
Cuspid with Distal Hook	0.46mm (.018")	L	297-418CV
Torque 0°		R	297-419CV
Tip +8°	0.56mm (.022")	L	297-420CV
		R	297-421CV
Universal Bicuspid with Hook	0.46mm (.018")	L	297-284CV
Torque -7°		R	297-285CV
Tip 0°	0.56mm (.022")	L	297-286CV
		R	297-287CV

Order individual brackets by placing a "U" in front of the item number.
 Package of 5

MANDIBULAR BRACKETS

	Slot Size	Left/Right	Bondable
Anteriors	0.46mm (.018")	L/R	297-169CV
Torque -1°			
Tip 0°	0.56mm (.022")	L/R	297-160CV
Cuspid with Distal Hook	0.46mm (.018")	L	297-426CV
Torque -11°		R	297-427CV
Tip +6°	0.56mm (.022")	L	297-428CV
		R	297-429CV
1st Bicuspid with Hook	0.46mm (.018")	L	297-288CV
Torque -17°		R	297-289CV
Tip 0°	0.56mm (.022")	L	297-290CV
		R	297-291CV
2nd Bicuspid with Hook	0.46mm (.018")	L	297-292CV
Torque -22°		R	297-293CV
Tip 0°	0.56mm (.022")	L	297-294CV
		R	297-295CV

Order individual brackets by placing a "U" in front of the item number.
 Package of 5

ClearVu® Mini Ceramic Brackets

MBT* McLaughlin, Bennett, Trevisi System

PATIENT KITS

	With Hooks	0.46mm (.018") Slot	0.56mm (.022") Slot
Kits contain all ceramic brackets			
Upper 3-3	Cuspid	297-380CV	297-381CV
Upper 5-5	Cuspid/Bicuspid	297-408CV	297-409CV
Upper & Lower 3-3	Cuspid	297-382CV	297-383CV
Upper & Lower 5-5	Cuspid/Bicuspid	297-410CV	297-411CV

Kits combine ceramic with metal* brackets

Up. 5-5 Ceramic with Lo. 5-5 Metal	Cuspid/Bicuspid	297-031CV	297-036CV
------------------------------------	-----------------	------------------	------------------

Kits contain maxillary first and second bicuspid with .027" in/out compensation.

*Nu-Edge Cobalt Chromium brackets.

297-495CV Bi-Dimensional Brackets with Hooks
.018" centrals and laterals with .022" cuspids and bicuspid
Kit contains 20 brackets.



Patented

MAXILLARY BRACKETS

	Slot Size	Left/Right	Bondable
Central	0.46mm	L	297-310CV
Torque +17°	(.018")	R	297-311CV
Tip +4°	0.56mm	L	297-312CV
	(.022")	R	297-313CV
Lateral	0.46mm	L	297-314CV
Torque +10°	(.018")	R	297-315CV
Tip +8°	0.56mm	L	297-316CV
	(.022")	R	297-317CV
Cuspid with Distal Hook	0.46mm	L	297-418CV
Torque 0°	(.018")	R	297-419CV
Tip +8°	0.56mm	L	297-420CV
	(.022")	R	297-421CV
Universal Bicuspid with Hook	0.46mm	L	297-285CV
Torque -7°	(.018")	R	297-284CV
Tip 0°	0.56mm	L	297-287CV
In/Out Comp. .027"	(.022")	R	297-286CV

Order individual brackets by placing a "U" in front of the item number.
Package of 5

MANDIBULAR BRACKETS

	Slot Size	Left/Right	Bondable
Anteriors Narrow	0.46mm	L/R	297-354CV
Torque -6°	(.018")		
Tip 0°	0.56mm	L/R	297-355CV
	(.022")		
Cuspid with Distal Hook	0.46mm	L	297-446CV
Torque 0°	(.018")	R	297-447CV
Tip +3°	0.56mm	L	297-448CV
	(.022")	R	297-449CV
1st Bicuspid with Hook	0.46mm	L	297-296CV
Torque -12°	(.018")	R	297-297CV
Tip 0°	0.56mm	L	297-298CV
	(.022")	R	297-299CV
2nd Bicuspid with Hook	0.46mm	L	297-288CV
Torque -17°	(.018")	R	297-289CV
Tip 0°	0.56mm	L	297-290CV
	(.022")	R	297-291CV

Order individual brackets by placing a "U" in front of the item number.
Package of 5


*MBT is a registered trademark of 3M. TP Orthodontics version of the McLaughlin, Bennett and Trevisi prescription is not intended to be a duplication of, nor an endorsement by, these doctors and such endorsement is not implied by the use of this term.

ClearVu® Mini Ceramic Brackets with Pre-Applied Adhesive

Compatible with the Roth system: not endorsed by Dr. Roth

PATIENT KITS

	With Hooks	0.46mm (.018") Slot	0.56mm (.022") Slot
Kits contain all ceramic brackets			
Upper 3-3	Cuspid	297-100CVX	297-101CVX
Upper 5-5	Cuspid/Bicuspid	297-404CVX	297-405CVX
Upper & Lower 3-3	Cuspid	297-102CVX	297-103CVX
Upper & Lower 5-5	Cuspid/Bicuspid	297-406CVX	297-407CVX



Patented

297-490CVX Bi-Dimensional Brackets with Hooks

.018" centrals and laterals with .022" cuspids and bicuspids
Kit contains 20 brackets.

MAXILLARY BRACKETS

	Slot Size	Left/Right	Bondable
Central	0.46mm	L	297-138CVX
Torque +14°	(.018")	R	297-139CVX
Tip +5°	0.56mm	L	297-150CVX
	(.022")	R	297-151CVX
Lateral	0.46mm	L	297-148CVX
Torque +7°	(.018")	R	297-149CVX
Tip +8°	0.56mm	L	297-152CVX
	(.022")	R	297-153CVX
Cuspid with Distal Hook	0.46mm	L	297-418CVX
Torque 0°	(.018")	R	297-419CVX
Tip +8°	0.56mm	L	297-420CVX
	(.022")	R	297-421CVX
Bicuspid with Hook	0.46mm	L	297-284CVX
Torque -7°	(.018")	R	297-285CVX
Tip 0°	0.56mm	L	297-286CVX
	(.022")	R	297-287CVX

Order individual brackets by placing a "U" in front of the item number.
Package of 5

MANDIBULAR BRACKETS

	Slot Size	Left/Right	Bondable
Anteriors	0.46mm	L/R	297-169CVX
Torque -1°	(.018")		
Tip 0°	0.56mm	L/R	297-160CVX
	(.022")		
Cuspid with Distal Hook	0.46mm	L	297-426CVX
Torque -11°	(.018")	R	297-427CVX
Tip +6°	0.56mm	L	297-428CVX
	(.022")	R	297-429CVX
1st Bicuspid with Hook	0.46mm	L	297-288CVX
Torque -17°	(.018")	R	297-289CVX
Tip 0°	0.56mm	L	297-290CVX
	(.022")	R	297-291CVX
2nd Bicuspid with Hook	0.46mm	L	297-292CVX
Torque -22°	(.018")	R	297-293CVX
Tip 0°	0.56mm	L	297-294CVX
	(.022")	R	297-295CVX

Order individual brackets by placing a "U" in front of the item number.
Package of 5


ClearVu® Mini Ceramic Brackets with Pre-Applied Adhesive

MBT* McLaughlin, Bennett, Trevisi System

PATIENT KITS

	With Hooks	0.46mm (.018") Slot	0.56mm (.022") Slot
Kits contain all ceramic brackets			
Upper 3-3	Cuspid	297-380CVX	297-381CVX
Upper 5-5	Cuspid/Bicuspid	297-408CVX	297-409CVX
Upper & Lower 3-3	Cuspid	297-382CVX	297-383CVX
Upper & Lower 5-5	Cuspid/Bicuspid	297-410CVX	297-411CVX

Kits contain maxillary first and second bicuspids with .027" in/out compensation.



Patented

297-495CVX Bi-Dimensional Brackets with Hooks

.018" centrals and laterals with .022" cuspids and bicuspids

Kit contains 20 brackets.

MAXILLARY BRACKETS

	Slot Size	Left/Right	Bondable
Central	0.46mm	L	297-310CVX
Torque +17°	(.018")	R	297-311CVX
Tip +4°	0.56mm	L	297-312CVX
	(.022")	R	297-313CVX
Lateral	0.46mm	L	297-314CVX
Torque +10°	(.018")	R	297-315CVX
Tip +8°	0.56mm	L	297-316CVX
	(.022")	R	297-317CVX
Cuspid with Distal Hook	0.46mm	L	297-418CVX
Torque 0°	(.018")	R	297-419CVX
Tip +8°	0.56mm	L	297-420CVX
	(.022")	R	297-421CVX
Bicuspid with Hook	0.46mm	L	297-285CVX
Torque -7°	(.018")	R	297-284CVX
Tip 0°	0.56mm	L	297-287CVX
In/Out Comp. .027"	(.022")	R	297-286CVX

Package of 5

MANDIBULAR BRACKETS

	Slot Size	Left/Right	Bondable
Anteriors Narrow	0.46mm	L/R	297-354CVX
Torque -6°	(.018")		
Tip 0°	0.56mm	L/R	297-355CVX
	(.022")		
Cuspid with Distal Hook	0.46mm	L	297-446CVX
Torque 0°	(.018")	R	297-447CVX
Tip +3°	0.56mm	L	297-448CVX
	(.022")	R	297-449CVX
1st Bicuspid with Hook	0.46mm	L	297-296CVX
Torque -12°	(.018")	R	297-297CVX
Tip 0°	0.56mm	L	297-298CVX
	(.022")	R	297-299CVX
2nd Bicuspid with Hook	0.46mm	L	297-288CVX
Torque -17°	(.018")	R	297-289CVX
Tip 0°	0.56mm	L	297-290CVX
In/Out Comp. .047"	(.022")	R	297-291CVX

Package of 5


*MBT is a registered trademark of 3M. TP Orthodontics version of the McLaughlin, Bennett and Trevisi prescription is not intended to be a duplication of, nor an endorsement by, these doctors and such endorsement is not implied by the use of this term.

CLEARVU®


Mini Ceramic Brackets



BRACKET ID CHART - ClearVu Roth



BRACKET ID CHART - ClearVu MBT



InVu® Mini Ceramic Brackets

Compatible with the Roth system: not endorsed by Dr. Roth

PATIENT KITS

	With Hooks	0.46mm (.018") Slot	0.56mm (.022") Slot
Kits contain all ceramic brackets			
Upper 3-3	No Hooks	297-106	297-110
Upper 5-5	Cuspid	297-060	297-065
Upper 5-5	Cuspid/Bicuspid	297-404	297-405
Upper 5-5	No Hooks	297-107	297-111
Upper & Lower 3-3	No Hooks	297-108	297-112
Upper & Lower 5-5	Cuspid	297-070	297-075
Upper & Lower 5-5	Cuspid/Bicuspid	297-406	297-407
Upper & Lower 5-5	No Hooks	297-109	297-113

Kits contain maxillary first and second bicuspids with .027" in/out compensation.



Patented

MAXILLARY BRACKETS

	Slot Size	Left/Right	Bondable
Central	0.46mm	L	297-138
Torque +14°	(.018")	R	297-139
Tip +5°	0.56mm	L	297-150
	(.022")	R	297-151
Lateral	0.46mm	L	297-148
Torque +7°	(.018")	R	297-149
Tip +8°	0.56mm	L	297-152
	(.022")	R	297-153
Cuspid without Hook	0.46mm	L	297-158
Torque 0°	(.018")	R	297-159
Tip +8°	0.56mm	L	297-154
	(.022")	R	297-155
Cuspid with Distal Hook	0.46mm	L	297-418
Torque 0°	(.018")	R	297-419
Tip +8°	0.56mm	L	297-420
	(.022")	R	297-421
Universal Bicuspid w/o Hook	0.46mm	L/R	297-226M
Torque -7°	(.018")		
Tip 0°	0.56mm	L/R	297-228M
	(.022")		
Universal Bicuspid w/ Hook	0.46mm	L	297-284M
Torque -7°	(.018")	R	297-285M
Tip 0°	0.56mm	L	297-286M
	(.022")	R	297-287M

Order individual brackets by placing a "U" in front of the item number.
Package of 5

MANDIBULAR BRACKETS

	Slot Size	Left/Right	Bondable
Anteriors	0.46mm	L/R	297-169
Torque -1°	(.018")		
Tip 0°	0.56mm	L/R	297-160
	(.022")		
Cuspid without Hook	0.46mm	L	297-176
Torque -11°	(.018")	R	297-177
Tip +6°	0.56mm	L	297-174
	(.022")	R	297-175
Cuspid with Distal Hook	0.46mm	L	297-426
Torque -11°	(.018")	R	297-427
Tip +6°	0.56mm	L	297-428
	(.022")	R	297-429
1st Bicuspid without Hook	0.46mm	L/R	297-182M
Torque -17°	(.018")		
Tip 0°	0.56mm	L/R	297-184M
	(.022")		
1st Bicuspid with Hook	0.46mm	L	297-288M
Torque -17°	(.018")	R	297-289M
Tip 0°	0.56mm	L	297-290M
	(.022")	R	297-291M
2nd Bicuspid without Hook	0.46mm	L/R	297-186M
Torque -22°	(.018")		
Tip 0°	0.56mm	L/R	297-188M
	(.022")		
2nd Bicuspid with Hook	0.46mm	L	297-292M
Torque -22°	(.018")	R	297-293M
Tip 0°	0.56mm	L	297-294M
	(.022")	R	297-295M

Order individual brackets by placing a "U" in front of the item number.
Package of 5

InVu® Ceramic Brackets

Compatible with the Roth system: not endorsed by Dr. Roth

PATIENT KITS

	With Hooks	0.46mm (.018") Slot	0.56mm (.022") Slot
Kits contain all ceramic brackets			
Upper 3-3	Cuspid	297-120	297-121
Upper 5-5	Cuspid	297-124	297-125
Upper 5-5	Cuspid/Bicuspid	297-126	297-127
Upper & Lower 3-3	Cuspid	297-132	297-133
Upper & Lower 5-5	Cuspid	297-136	297-137
Upper & Lower 5-5	Cuspid/Bicuspid	297-128	297-129
Kits combine ceramic with metal* brackets			
Up. 3-3 Ceramic & Up. 4-5 Metal with Lo. 5-5 Metal	Cuspid	297-040	297-045
Up. 4-4 Ceramic & Up. 5 Metal with Lo. 3-3 Ceramic & Lo. 4-5 Metal	Cuspid	297-056	297-057
Up. 5-5 Ceramic with Lo. 5-5 Metal	Cuspid	297-050	297-055

*Nu-Edge Cobalt Chromium brackets.



Patented

MAXILLARY BRACKETS

	Slot Size	Left/Right	Bondable
Central	0.46mm	L	297-140
Torque +14°	(.018")	R	297-141
Tip +5°	0.56mm	L	297-142
	(.022")	R	297-143
Lateral	0.46mm	L	297-144
Torque +7°	(.018")	R	297-145
Tip +8°	0.56mm	L	297-146
	(.022")	R	297-147
Cuspid with Distal Hook	0.46mm	L	297-222
Torque 0°	(.018")	R	297-223
Tip +8°	0.56mm	L	297-224
	(.022")	R	297-225
Universal Bicuspid w/o Hook	0.46mm	L/R	297-226
Torque -7°	(.018")		
Tip 0°	0.56mm	L/R	297-228
	(.022")		
Universal Bicuspid with Hook	0.46mm	L	297-284
Torque -7°	(.018")	R	297-285
Tip 0°	0.56mm	L	297-286
	(.022")	R	297-287

Order individual brackets by placing a "U" in front of the item number.

Package of 5

MANDIBULAR BRACKETS

	Slot Size	Left/Right	Bondable
Anteriors	0.46mm	L/R	297-170
Torque -1°	(.018")		
Tip 0°	0.56mm	L/R	297-172
	(.022")		
Cuspid with Distal Hook	0.46mm	L	297-178
Torque -11°	(.018")	R	297-179
Tip +6°	0.56mm	L	297-180
	(.022")	R	297-181
1st Bicuspid without Hook	0.46mm	L/R	297-182
Torque -17°	(.018")		
Tip 0°	0.56mm	L/R	297-184
	(.022")		
1st Bicuspid with Hook	0.46mm	L	297-288
Torque -17°	(.018")	R	297-289
Tip 0°	0.56mm	L	297-290
	(.022")	R	297-291
2nd Bicuspid without Hook	0.46mm	L/R	297-186
Torque -22°	(.018")		
Tip 0°	0.56mm	L/R	297-188
	(.022")		
2nd Bicuspid with Hook	0.46mm	L	297-292
Torque -22°	(.018")	R	297-293
Tip 0°	0.56mm	L	297-294
	(.022")	R	297-295

Order individual brackets by placing a "U" in front of the item number.

Package of 5

InVu® Mini Ceramic Brackets

MBT* McLaughlin, Bennett, Trevisi System

PATIENT KITS

	With Hooks	0.46mm (.018") Slot	0.56mm (.022") Slot
Kits contain all ceramic brackets			
Upper 3-3	No Hooks	297-391	297-393
Upper 5-5	Cuspid	297-394	297-396
Upper 5-5	Cuspid/Bicuspid	297-408	297-409
Upper 5-5	No Hooks	297-395	297-397
Upper & Lower 3-3	No Hooks	297-303	297-304
Upper & Lower 5-5	Cuspid	297-306	297-308
Upper & Lower 5-5	Cuspid/Bicuspid	297-410	297-411
Upper & Lower 5-5	No Hooks	297-307	297-309

Kits contain maxillary first and second bicuspids with .027" in/out compensation.



Patented

MAXILLARY BRACKETS

	Slot Size	Left/Right	Bondable
Central	0.46mm	L	297-310
Torque +17°	(.018")	R	297-311
Tip +4°	0.56mm	L	297-312
	(.022")	R	297-313
Lateral	0.46mm	L	297-314
Torque +10°	(.018")	R	297-315
Tip +8°	0.56mm	L	297-316
	(.022")	R	297-317
Cuspid without Hook	0.46mm	L	297-318
Torque -7°	(.018")	R	297-319
Tip +8°	0.56mm	L	297-320
	(.022")	R	297-321
Cuspid with Distal Hook	0.46mm	L	297-418
Torque 0°	(.018")	R	297-419
Tip +8°	0.56mm	L	297-420
	(.022")	R	297-421
Universal Bicuspid w/o Hook	0.46mm	L/R	297-226M
Torque -7°	(.018")		
Tip 0°	0.56mm	L/R	297-228M
	(.022")		
Universal Bicuspid w Hook	0.46mm	L	297-285M
Torque -7°	(.018")	R	297-284M
Tip 0°	0.56mm	L	297-287M
	(.022")	R	297-286M

Order individual brackets by placing a "U" in front of the item number.

Package of 5

MANDIBULAR BRACKETS

	Slot Size	Left/Right	Bondable
Anteriors	0.46mm	L/R	297-354
Torque -6°	(.018")		
Tip 0°	0.56mm	L/R	297-355
	(.022")		
Cuspid without Hook	0.46mm	L	297-356
Torque -6°	(.018")	R	297-357
Tip +3°	0.56mm	L	297-358
	(.022")	R	297-359
Cuspid with Distal Hook	0.46mm	L	297-446
Torque 0°	(.018")	R	297-447
Tip +3°	0.56mm	L	297-448
	(.022")	R	297-449
1st Bicuspid without Hook	0.46mm	L/R	297-264M
Torque -12°	(.018")		
Tip 0°	0.56mm	L/R	297-266M
	(.022")		
1st Bicuspid with Hook	0.46mm	L	297-296M
Torque -12°	(.018")	R	297-297M
Tip 0°	0.56mm	L	297-298M
	(.022")	R	297-299M
2nd Bicuspid without Hook	0.46mm	L/R	297-272M
Torque -17°	(.018")		
Tip 0°	0.56mm	L/R	297-274M
	(.022")		
2nd Bicuspid with Hook	0.46mm	L	297-288M
Torque -17°	(.018")	R	297-289M
Tip 0°	0.56mm	L	297-290M
	(.022")	R	297-291M

Order individual brackets by placing a "U" in front of the item number.

Package of 5

*MBT is a registered trademark of 3M. TP Orthodontics version of the McLaughlin, Bennett and Trevisi prescription is not intended to be a duplication of, nor an endorsement by, these doctors and such endorsement is not implied by the use of this term.

InVu® Ceramic Brackets


MBT* McLaughlin, Bennett, Trevisi System

PATIENT KITS

	With Hooks	0.46mm (.018") Slot	0.56mm (.022") Slot
Kits contain all ceramic brackets			
Upper 3-3	Cuspid	297-190	297-192
	No Hooks	297-191	297-193
Upper 5-5	Cuspid	297-194	297-196
	Cuspid/Bicuspid	297-400	297-401
Upper & Lower 3-3	No Hooks	297-195	297-197
	Cuspid	297-202	297-204
Upper & Lower 5-5	No Hooks	297-203	297-205
	Cuspid	297-206	297-208
	Cuspid/Bicuspid	297-402	297-403
	No Hooks	297-207	297-209
Kits combine ceramic with metal* brackets			
Up. 3-3 Ceramic & Up. 4-5 Metal with Lo. 5-5 Metal	Cuspid	297-020	297-025
Up. 4-4 Ceramic & Up. 5-5 Metal with Lo. 3-3 Ceramic & Lo. 4-5 Metal	Cuspid	297-036	297-037
Up. 4-4 Ceramic & Up. 5-5 Metal with Lo. 3-3 Ceramic & Lo. 4-5 Metal	No Hooks	297-038	297-039
Up. 5-5 Ceramic with Lo. 5-5 Metal	Cuspid	297-030	297-035

Kits contain maxillary first and second bicuspids with .027" in/out compensation.

*Nu-Edge Cobalt Chromium brackets.



Patented

MAXILLARY BRACKETS

	Slot Size	Left/ Right	Bondable
Central	0.46mm	L	297-210
Torque +17°	(.018")	R	297-211
Tip +4°	0.56mm	L	297-212
	(.022")	R	297-213
Lateral	0.46mm	L	297-214
Torque +10°	(.018")	R	297-215
Tip +8°	0.56mm	L	297-216
	(.022")	R	297-217
Cuspid without Hook	0.46mm	L	297-218
Torque -7°	(.018")	R	297-219
Tip +8°	0.56mm	L	297-220
	(.022")	R	297-221
Cuspid with Distal Hook	0.46mm	L	297-222
Torque 0°	(.018")	R	297-223
Tip +8°	0.56mm	L	297-224
	(.022")	R	297-225
Universal Bicuspid without Hook	0.46mm	L/R	297-226
Torque -7°	(.018")		
Tip 0°	0.56mm	L/R	297-228
In/Out Comp. .027"	(.022")		
Universal Bicuspid with Hook	0.46mm	L	297-285
Torque -7°	(.018")	R	297-284
Tip 0°	0.56mm	L	297-287
In/Out Comp. .027"	(.022")	R	297-286

Order individual brackets by placing a "U" in front of the item number.
Package of 5

MANDIBULAR BRACKETS

	Slot Size	Left/ Right	Bondable
Anteriors Narrow	0.46mm	L/R	297-254
Torque -6°	(.018")		
Tip 0°	0.56mm	L/R	297-255
	(.022")		
Anteriors Wide	0.46mm	L/R	297-236
Torque -6°	(.018")		
Tip 0°	0.56mm	L/R	297-237
	(.022")		
Cuspid without Hook	0.46mm	L	297-256
Torque -6°	(.018")	R	297-257
Tip +3°	0.56mm	L	297-258
	(.022")	R	297-259
Cuspid with Distal Hook	0.46mm	L	297-260
Torque 0°	(.018")	R	297-261
Tip +3°	0.56mm	L	297-262
	(.022")	R	297-263
1st Bicuspid without Hook	0.46mm	L/R	297-264
Torque -12°	(.018")		
Tip 0°	0.56mm	L/R	297-266
	(.022")		
1st Bicuspid with Hook	0.46mm	L	297-296
Torque -12°	(.018")	R	297-297
Tip 0°	0.56mm	L	297-298
	(.022")	R	297-299
2nd Bicuspid without Hook	0.46mm	L/R	297-272
Torque -17°	(.018")		
Tip 0°	0.56mm	L/R	297-274
	(.022")		
2nd Bicuspid with Hook	0.46mm	L	297-288
Torque -17°	(.018")	R	297-289
Tip 0°	0.56mm	L	297-290
	(.022")	R	297-291


Order individual brackets by placing a "U" in front of the item number.
Package of 5

InVu® Ceramic Brackets

Standard Edgewise System

PATIENT KITS

	0.46mm (.018") Slot	0.56mm (.022") Slot
Kits contain all ceramic brackets without cuspid hooks		
Upper 3-3	297-960	297-961
Upper 5-5	297-962	297-963
Upper & Lower 5-5	297-964	297-965



Patented

MAXILLARY BRACKETS

	Slot Size	Left/ Right	Bondable
Central	0.46mm	L/R	297-905
Torque +0°	(.018")		
Tip +0°	0.56mm	L/R	297-906
	(.022")		
Lateral	0.46mm	L/R	297-911
Torque +0°	(.018")		
Tip +0°	0.56mm	L/R	297-912
	(.022")		
Cuspid	0.46mm	L/R	297-929
Torque +0°	(.018")		
Tip +0°	0.56mm	L/R	297-930
	(.022")		
Bicuspid without Hook	0.46mm	L/R	297-939
Torque +0°	(.018")		
Tip +0°	0.56mm	L/R	297-940
	(.022")		

Order individual brackets by placing a "U" in front of the item number.

Package of 5

MANDIBULAR BRACKETS

	Slot Size	Left/ Right	Bondable
Anteriors	0.46mm	L/R	297-933
Torque +0°	(.018")		
Tip +0°	0.56mm	L/R	297-934
	(.022")		
Cuspid	0.46mm	L/R	297-929
Torque +0°	(.018")		
Tip +0°	0.56mm	L/R	297-930
	(.022")		
Bicuspid without Hook	0.46mm	L/R	297-941
Torque +0°	(.018")		
Tip +0°	0.56mm	L/R	297-942
	(.022")		

Order individual brackets by placing a "U" in front of the item number.


Package of 5

InVu® Mini Ceramic Brackets with Pre-Applied Adhesive

Compatible with the Roth system: not endorsed by Dr. Roth

PATIENT KITS

	With Hooks	0.46mm (.018") Slot	0.56mm (.022") Slot
Kits contain all ceramic brackets			
Upper 5-5	Cuspid	297-060RBX	297-065RBX
Upper 5-5	Cuspid/Bicuspid	297-404RBX	297-405RBX
Upper 5-5	No Hooks	297-107RBX	297-111RBX
Upper & Lower 5-5	Cuspid	297-070RBX	297-075RBX
Upper & Lower 5-5	Cuspid/Bicuspid	297-406RBX	297-407RBX
Upper & Lower 5-5	No Hooks	297-109RBX	297-113RBX



Patented

MAXILLARY BRACKETS — MINI

	Slot Size	Left/Right	Bondable
Central	0.46mm	L	297-138RBX
Torque +14°	(.018")	R	297-139RBX
Tip +5°	0.56mm	L	297-150RBX
	(.022")	R	297-151RBX
Lateral	0.46mm	L	297-148RBX
Torque +7°	(.018")	R	297-149RBX
Tip +8°	0.56mm	L	297-152RBX
	(.022")	R	297-153RBX
Cuspid without Hook	0.46mm	L	297-158RBX
Torque +0°	(.018")	R	297-159RBX
Tip +8°	0.56mm	L	297-154RBX
	(.022")	R	297-155RBX
Cuspid with Distal Hook	0.46mm	L	297-418RBX
Torque +0°	(.018")	R	297-419RBX
Tip +8°	0.56mm	L	297-420RBX
	(.022")	R	297-421RBX
Universal Bicuspid w/o Hook	0.46mm	L/R	297-226MRX
Torque -7°	(.018")		
Tip 0°	0.56mm	L/R	297-228MRX
	(.022")		
Bicuspid with Hook	0.46mm	L	297-284MRX
Torque -7°	(.018")	R	297-285MRX
Tip 0°	0.56mm	L	297-286MRX
	(.022")	R	297-287MRX

Package of 5

MANDIBULAR BRACKETS — MINI

	Slot Size	Left/Right	Bondable
Anteriors	0.46mm	L/R	297-169RBX
Torque -1°	(.018")		
Tip 0°	0.56mm	L/R	297-160RBX
	(.022")		
Cuspid without Hook	0.46mm	L	297-176RBX
Torque -11°	(.018")	R	297-177RBX
Tip +6°	0.56mm	L	297-174RBX
	(.022")	R	297-175RBX
Cuspid with Distal Hook	0.46mm	L	297-426RBX
Torque -11°	(.018")	R	297-427RBX
Tip +6°	0.56mm	L	297-428RBX
	(.022")	R	297-429RBX
1st Bicuspid without Hook	0.46mm	L/R	297-182MRX
Torque -17°	(.018")		
Tip 0°	0.56mm	L/R	297-184MRX
	(.022")		
1st Bicuspid with Hook	0.46mm	L	297-288MRX
Torque -17°	(.018")	R	297-289MRX
Tip 0°	0.56mm	L	297-290MRX
	(.022")	R	297-291MRX
2nd Bicuspid without Hook	0.46mm	L/R	297-186MRX
Torque -22°	(.018")		
Tip 0°	0.56mm	L/R	297-188MRX
	(.022")		
2nd Bicuspid with Hook	0.46mm	L	297-292MRX
Torque -22°	(.018")	R	297-293MRX
Tip 0°	0.56mm	L	297-294MRX
	(.022")	R	297-295MRX


Package of 5

InVu® Ceramic Brackets with Pre-Applied Adhesive

Compatible with the Roth system: not endorsed by Dr. Roth

PATIENT KITS

	With Hooks	0.46mm (.018") Slot	0.56mm (.022") Slot
Kits contain all ceramic brackets			
Upper 5-5	Cuspid	297-124RBX	297-125RBX
Upper 5-5	Cuspid/Bicuspid	297-126RBX	297-127RBX
Upper & Lower 5-5	Cuspid	297-136RBX	297-137RBX
Upper & Lower 5-5	Cuspid/Bicuspid	297-128RBX	297-129RBX



Patented

MAXILLARY BRACKETS

	Slot Size	Left/ Right	Bondable
Central	0.46mm	L	297-140RBX
Torque +14°	(.018")	R	297-141RBX
Tip +5°	0.56mm	L	297-142RBX
	(.022")	R	297-143RBX
Lateral	0.46mm	L	297-144RBX
Torque +7°	(.018")	R	297-145RBX
Tip +8°	0.56mm	L	297-146RBX
	(.022")	R	297-147RBX
Cuspid with Distal Hook	0.46mm	L	297-222RBX
Torque 0°	(.018")	R	297-223RBX
Tip +8°	0.56mm	L	297-224RBX
	(.022")	R	297-225RBX
Universal Bicuspid w/o Hook	0.46mm	L/R	297-226RBX
Torque -7°	(.018")		
Tip 0°	0.56mm	L/R	297-228RBX
	(.022")		
Bicuspid with Hook	0.46mm	L	297-284RBX
Torque -7°	(.018")	R	297-285RBX
Tip 0°	0.56mm	L	297-286RBX
	(.022")	R	297-287RBX

Package of 5

MANDIBULAR BRACKETS

	Slot Size	Left/ Right	Bondable
Anteriors	0.46mm	L/R	297-170RBX
Torque -1°	(.018")		
Tip 0°	0.56mm	L/R	297-172RBX
	(.022")		
Cuspid with Distal Hook	0.46mm	L	297-178RBX
Torque -11°	(.018")	R	297-179RBX
Tip +6°	0.56mm	L	297-180RBX
	(.022")	R	297-181RBX
1st Bicuspid without Hook	0.46mm	L/R	297-182RBX
Torque -17°	(.018")		
Tip 0°	0.56mm	L/R	297-184RBX
	(.022")		
1st Bicuspid with Hook	0.46mm	L	297-288RBX
Torque -17°	(.018")	R	297-289RBX
Tip 0°	0.56mm	L	297-290RBX
	(.022")	R	297-291RBX
2nd Bicuspid without Hook	0.46mm	L/R	297-186RBX
Torque -22°	(.018")		
Tip 0°	0.56mm	L/R	297-188RBX
	(.022")		
2nd Bicuspid with Hook	0.46mm	L	297-292RBX
Torque -22°	(.018")	R	297-293RBX
Tip 0°	0.56mm	L	297-294RBX
	(.022")	R	297-295RBX


Package of 5

InVu® Mini Ceramic Brackets with Pre-Applied Adhesive

MBT* McLaughlin, Bennett, Trevisi System

PATIENT KITS

	With Hooks	0.46mm (.018") Slot	0.56mm (.022") Slot
Kits contain all ceramic brackets			
Upper 3-3	No Hooks	297-391RBX	297-393RBX
Upper 5-5	Cuspid	297-394RBX	297-396RBX
Upper 5-5	Cuspid/Bicuspid	297-408RBX	297-409RBX
Upper 5-5	No Hooks	297-395RBX	297-397RBX
Upper & Lower 3-3	No Hooks	297-303RBX	297-304RBX
Upper & Lower 5-5	Cuspid	297-306RBX	297-308RBX
Upper & Lower 5-5	Cuspid/Bicuspid	297-410RBX	297-411RBX
Upper & Lower 5-5	No Hooks	297-307RBX	297-309RBX



Patented

MAXILLARY BRACKETS

	Slot Size	Left/ Right	Bondable
Central	0.46mm	L	297-310RBX
Torque +17°	(.018")	R	297-311RBX
Tip +4°	0.56mm	L	297-312RBX
	(.022")	R	297-313RBX
Lateral	0.46mm	L	297-314RBX
Torque +10°	(.018")	R	297-315RBX
Tip +8°	0.56mm	L	297-316RBX
	(.022")	R	297-317RBX
Cuspid without Hook	0.46mm	L	297-318RBX
Torque -7°	(.018")	R	297-319RBX
Tip +8°	0.56mm	L	297-320RBX
	(.022")	R	297-321RBX
Cuspid with Distal Hook	0.46mm	L	297-418RBX
Torque +0°	(.018")	R	297-419RBX
Tip +8°	0.56mm	L	297-420RBX
	(.022")	R	297-421RBX
Universal Bicuspid w/o Hook	0.46mm	L/R	297-226MRX
Torque -7°	(.018")		
Tip 0°	0.56mm	L/R	297-228MRX
	(.022")		
Bicuspid with Hook	0.46mm	L	297-285MRX
Torque -7°	(.018")	R	297-284MRX
Tip 0°	0.56mm	L	297-287MRX
	(.022")	R	297-286MRX

Package of 5

MANDIBULAR BRACKETS

	Slot Size	Left/ Right	Bondable
Anteriors	0.46mm	L/R	297-354RBX
Torque -6°	(.018")		
Tip 0°	0.56mm	L/R	297-355RBX
	(.022")		
Cuspid without Hook	0.46mm	L	297-356RBX
Torque -6°	(.018")	R	297-357RBX
Tip +3°	0.56mm	L	297-358RBX
	(.022")	R	297-359RBX
Cuspid with Distal Hook	0.46mm	L	297-446RBX
Torque 0°	(.018")	R	297-447RBX
Tip +3°	0.56mm	L	297-448RBX
	(.022")	R	297-449RBX
1st Bicuspid without Hook	0.46mm	L/R	297-264MRX
Torque -12°	(.018")		
Tip 0°	0.56mm	L/R	297-266MRX
	(.022")		
1st Bicuspid with Hook	0.46mm	L	297-296MRX
Torque -12°	(.018")	R	297-297MRX
Tip 0°	0.56mm	L	297-298MRX
	(.022")	R	297-299MRX
2nd Bicuspid without Hook	0.46mm	L/R	297-272MRX
Torque -17°	(.018")		
Tip 0°	0.56mm	L/R	297-274MRX
	(.022")		
2nd Bicuspid with Hook	0.46mm	L	297-288MRX
Torque -17°	(.018")	R	297-289MRX
Tip 0°	0.56mm	L	297-290MRX
	(.022")	R	297-291MRX

Package of 5

*MBT is a registered trademark of 3M. TP Orthodontics version of the McLaughlin, Bennett and Trevisi prescription is not intended to be a duplication of, nor an endorsement by, these doctors and such endorsement is not implied by the use of this term.


InVu® Ceramic Brackets with Pre-Applied Adhesive

MBT* McLaughlin, Bennett, Trevisi System

PATIENT KITS

	With Hooks	0.46mm (.018") Slot	0.56mm (.022") Slot
Kits contain all ceramic brackets			
Upper 5-5	Cuspid	297-194RBX	297-196RBX
	Cuspid/Bicuspid	297-400RBX	297-401RBX
Upper & Lower 5-5	Cuspid	297-206RBX	297-208RBX
	Cuspid/Bicuspid	297-402RBX	297-403RBX
	No Hooks	297-207RBX	297-209RBX

Kits contain maxillary first and second bicuspids with .027" in/out compensation.



Patented

MAXILLARY BRACKETS

	Slot Size	Left/ Right	Bondable
Central	0.46mm	L	297-210RBX
Torque +17°	(.018")	R	297-211RBX
Tip +4°	0.56mm	L	297-212RBX
	(.022")	R	297-213RBX
Lateral	0.46mm	L	297-214RBX
Torque +10°	(.018")	R	297-215RBX
Tip +8°	0.56mm	L	297-216RBX
	(.022")	R	297-217RBX
Cuspid without Hook	0.46mm	L	297-218RBX
Torque -7°	(.018")	R	297-219RBX
Tip +8°	0.56mm	L	297-220RBX
	(.022")	R	297-221RBX
Cuspid with Distal Hook	0.46mm	L	297-222RBX
Torque 0°	(.018")	R	297-223RBX
Tip +8°	0.56mm	L	297-224RBX
	(.022")	R	297-225RBX
Universal Bicuspid w/o Hook	0.46mm	L/R	297-226RBX
Torque -7°	(.018")		
Tip 0°	0.56mm	L/R	297-228RBX
In/Out Comp. .027"	(.022")		
Bicuspid with Hook	0.46mm	L	297-285RBX
Torque -7°	(.018")	R	297-284RBX
Tip 0°	0.56mm	L	297-287RBX
In/Out Comp. .027"	(.022")	R	297-286RBX

Package of 5

MANDIBULAR BRACKETS

	Slot Size	Left/ Right	Bondable
Anteriors Narrow	0.46mm	L/R	297-254RBX
Torque -6°	(.018")		
Tip 0°	0.56mm	L/R	297-255RBX
	(.022")		
Anteriors Wide	0.46mm	L/R	297-236RBX
Torque -6°	(.018")		
Tip 0°	0.56mm	L/R	297-237RBX
	(.022")		
Cuspid without Hook	0.46mm	L	297-256RBX
Torque -6°	(.018")	R	297-257RBX
Tip +3°	0.56mm	L	297-258RBX
	(.022")	R	297-259RBX
Cuspid with Distal Hook	0.46mm	L	297-260RBX
Torque 0°	(.018")	R	297-261RBX
Tip +3°	0.56mm	L	297-262RBX
	(.022")	R	297-263RBX
1st Bicuspid without Hook	0.46mm	L/R	297-264RBX
Torque -12°	(.018")		
Tip 0°	0.56mm	L/R	297-266RBX
	(.022")		
1st Bicuspid with Hook	0.46mm	L	297-296RBX
Torque -12°	(.018")	R	297-297RBX
Tip 0°	0.56mm	L	297-298RBX
	(.022")	R	297-299RBX
2nd Bicuspid without Hook	0.46mm	L/R	297-272RBX
Torque -17°	(.018")		
Tip 0°	0.56mm	L/R	297-274RBX
	(.022")		
2nd Bicuspid with Hook	0.46mm	L	297-288RBX
Torque -17°	(.018")	R	297-289RBX
Tip 0°	0.56mm	L	297-290RBX
	(.022")	R	297-291RBX

Package of 5


*MBT is a registered trademark of 3M. TP Orthodontics version of the McLaughlin, Bennett and Trevisi prescription is not intended to be a duplication of, nor an endorsement by, these doctors and such endorsement is not implied by the use of this term.




Ceramic Brackets




BRACKET ID CHART - InVu Roth



BRACKET ID CHART - InVu MBT



BRACKET ID CHART - InVu Standard Edgewise



Click-It® Self-Ligating Aesthetic Brackets

MBT* McLaughlin, Bennett, Trevisi System

PATIENT KITS

	With Cuspid Hooks	0.56mm (.022") Slot
Upper 5-5	Yes	299-000
	No	299-001
Upper & Lower 5-5	Yes	299-002
	No	299-003



Patented

MAXILLARY BRACKETS

	Slot Size	Left/ Right	Bondable
Central	0.56mm	L	299-012
Torque +17°	(.022")	R	299-013
Tip +4°			
Lateral	0.56mm	L	299-016
Torque +10°	(.022")	R	299-017
Tip +8°			
Cuspid without Hook	0.56mm	L	299-020
Torque 0°	(.022")	R	299-021
Tip +8°			
Cuspid with Distal Hook	0.56mm	L	299-024
Torque 0°	(.022")	R	299-025
Tip +8°			
Universal Bicuspid without Hook	0.56mm L/R		299-028
Torque -7°	(.022")		
Tip 0°			

Package of 5

MANDIBULAR BRACKETS

	Slot Size	Left/ Right	Bondable
Anteriors Narrow	0.56mm	L/R	299-055
Torque -6°	(.022")		
Tip 0°			
Cuspid without Hook	0.56mm	L	299-058
Torque 0°	(.022")	R	299-059
Tip +3°			
Cuspid with Distal Hook	0.56mm	L	299-062
Torque 0°	(.022")	R	299-063
Tip +3°			
1st Bicuspid without Hook	0.56mm	L	299-066
Torque -12°	(.022")	R	299-067
Tip 0°			
2nd Bicuspid without Hook	0.56mm	L	299-074
Torque -17°	(.022")	R	299-075
Tip 0°			

Package of 5

Forceps

Note: Forceps are required for use of the Click-It Self Ligating Bracket System. More information can be found on page 163.

100-198 Click-It Forceps Assembly

Each

100-098 Click-It Forceps Assembly, Angled, includes release tips

Each

100-199 Click-It Release Tips (5 Left, 5 Right)

Package of 10

100-200 Click-It Insertion Tips

Package of 5

100-099 Click-It Release Tips, Angled, 45° (5 Left, 5 Right)

Package of 10



*MBT is a registered trademark of 3M. TP Orthodontics version of the McLaughlin, Bennett and Trevisi prescription is not intended to be a duplication of, nor an endorsement by, these doctors and such endorsement is not implied by the use of this term.

CLICK-IT[®]

Aesthetic Self-Ligating Brackets



BRACKET ID CHART - Click-It MBT

RIGHT						LEFT					
2nd Bicuspid	1st Bicuspid	Cuspid	Lateral	Central		Central	Lateral	Cuspid	1st Bicuspid	2nd Bicuspid	
2nd Bicuspid	1st Bicuspid	Cuspid	Anteriors			Anteriors		Cuspid	1st Bicuspid	2nd Bicuspid	

Tip-Edge PLUS[®] Ceramic Brackets

PATIENT KITS

	Slot Size	Non-Oriented	Oriented Incisal Jig
Kits contain all ceramic brackets			
Upper 3-3	0.56mm	297-675	297-598
Upper & Lower 3-3	(.022")		297-599
Upper 5-5 Non-extraction		297-677	297-600
Upper & Lower 5-5 Non-extraction		297-678	297-602
Lower 3-3			297-604
Lower 5-5 Non-extraction			297-605

	Slot Size	Non-Oriented	Oriented Incisal Jig
Kits combine ceramic with metal brackets			
Up. 3-3 Ceramic & Up. 4-5 Metal with Lo. 5-5 Metal Non-extraction	0.56mm (.022")	297-608	297-610
Up. 5-5 Ceramic with Lo. 5-5 Metal Non-extraction	0.56mm (.022")	297-612	297-614

Non-extraction kits contain 20 bondable brackets, extraction kits have 16 brackets.

MAXILLARY BRACKETS

	Slot Size	Left/Right	Non-Oriented	Oriented 3.5mm Blue	Oriented 4mm Red	Oriented 4.5mm Yellow	Oriented 5mm Green
Central							
Torque +12°	0.56mm	L	297-611N		297-639	297-655	
Tip +5°	(.022")	R	297-612N		297-640	297-656	
Lateral							
Torque +8°	0.56mm	L	297-613N	297-629	297-641		
Tip +9°	(.022")	R	297-614N	297-630	297-642		
Cuspid							
Torque -2°	0.56mm	L	297-615		297-643	297-659	297-671
Tip +11°	(.022")	R	297-616		297-644	297-660	297-672
Bicuspid							
Torque -7°	0.56mm	CCW	297-617	297-631	297-645	297-661	
Tip 0°	(.022")	CW	297-618	297-632	297-646	297-662	

Package of 5

MANDIBULAR BRACKETS

Anterior							
Torque -1°	0.56mm	L	297-619	297-633	297-647		
Tip 0°	(.022")	R	297-620	297-634	297-648		
Cuspid							
Torque -11°	0.56mm	L	297-623		297-651	297-667	297-673
Tip +5°	(.022")	R	297-624		297-652	297-668	297-674
Bicuspid							
Torque -20°	0.56mm	CCW	297-625	297-637	297-653	297-669	
Tip 0°	(.022")	CW	297-626	297-638	297-654	297-670	

Recommended jig heights: Maxillary central 4mm, lateral 3.5mm, cuspid 4.5mm, bicuspid 4mm.
Mandibular anteriors 3.5mm, cuspids 4mm, bicuspid 4mm.

Package of 5

Tip-Edge PLUS[®] Ceramic Brackets with Pre-Applied Adhesive

PATIENT KITS

	Slot Size	Non-Oriented	Oriented Incisal Jig
Kits contain all ceramic brackets			
Upper & Lower 5-5	0.56mm	RB297-678	RB297-602
Non-extraction	(.022")		



Patented

MAXILLARY BRACKETS

	Slot Size	Left/Right	Non-Oriented	Oriented 3.5mm Blue	Oriented 4mm Red	Oriented 4.5mm Yellow
Central						
Torque +12°	0.56mm	L	RB297-611N		RB297-639	
Tip +5°	(.022")	R	RB297-612N		RB297-640	
Lateral						
Torque +8°	0.56mm	L	RB297-613N	RB297-629		
Tip +9°	(.022")	R	RB297-614N	RB297-630		
Cuspid						
Torque -2°	0.56mm	L	RB297-615			RB297-659
Tip +11°	(.022")	R	RB297-616			RB297-660
Bicuspid						
Torque -7°	0.56mm	CCW	RB297-617		RB297-645	
Tip 0°	(.022")	CW	RB297-618		RB297-646	

Package of 5

MANDIBULAR BRACKETS

Anterior						
Torque -1°	0.56mm	L	RB297-619	RB297-633		
Tip 0°	(.022")	R	RB297-620	RB297-634		
Cuspid						
Torque -11°	0.56mm	L	RB297-623		RB297-651	
Tip +5°	(.022")	R	RB297-624		RB297-652	
Bicuspid						
Torque -20°	0.56mm	CCW	RB297-625		RB297-653	
Tip 0°	(.022")	CW	RB297-626		RB297-654	

Recommended jig heights: Maxillary central 4mm, lateral 3.5mm, cuspid 4.5mm, bicuspid 4mm.
Mandibular anteriors 3.5mm, cuspids 4mm, bicuspid 4mm.


Package of 5

Tip-Edge PLUS®


Ceramic Brackets



BRACKET ID CHART - Tip-Edge PLUS Non-Extraction



BRACKET ID CHART - Tip-Edge PLUS Extraction




Circles identify the upper brackets and triangles identify the lower brackets. Arrows on the bicuspid represent the direction crowns should tip during initial treatment – either clockwise or counterclockwise. The bonding pad shape also points toward the direction of desired crown tipping.


Positioning Jigs

Incisal jigs are designed to help with accurate incisogingival bracket placement and assist with sighting long axis crown angulations.

Jig measurements are from the center of the archwire slot to the incisal edge or cusp tips. Color coding identifies the various heights.



RECOMMENDED JIG HEIGHT



Jig height substitutions can be made for an additional charge.
Green 5 mm jigs are also available for cuspids.

Gingivally offset mandibular bicuspid brackets are available with brackets positioned gingivally on bonding pads to permit proper bonding height on partially erupted teeth

Notes:

[illegible]

Visit us online at **tportho.com**