

Nu-Edge NX[®]

Cobalt Chromium Pre-Coated Brackets

Pre-coated for convenience.

Biocompatible for nickel sensitive patients.

Nu-Edge NX brackets combine the performance of a cobalt chromium bracket with the convenience of pre-applied adhesive, providing efficiency and results you can count on.

- Pre-coated with eXact orthodontic adhesive for easy clean-up
- Exceptionally strong CoCr helps prevent archwire slot collapse
- Polymer interface between the bracket and adhesive for enhanced bond strength

MBT PATIENT KITS

Slot Size	With Hooks	
0.46mm (.018")	Cuspid/Bicuspid	293-381NX
0.56mm (.022")	Cuspid/Bicuspid	293-384NX

Kits contain 20 brackets.

ROTH PATIENT KITS

Slot Size	With Hooks	
0.46mm (.018")	Cuspid/Bicuspid	293-206NX
0.56mm (.022")	Cuspid/Bicuspid	293-006NX

Kits contain 20 brackets.

Save Time with Pre-Applied Adhesive

Nu-Edge pre-coated brackets feature Readibase[®], with eXact[®] pre-applied adhesive for convenience and increased efficiency.

- Ready to place right out of the package
- Reduced risk of contamination
- Brackets set easily without sliding
- Minimal flash clean-up
- Reduction in bonding appointment time



Find out more at tportho.com
or call **800-348-8856**



TP Orthodontics, Inc.
tportho.com

Ultra Low Profile

Molar Tube


A closer fit for a stronger bond.

Positions easily and accurately.

Conforms to buccal groove.

Improves patient comfort.

Our Ultra Low Profile Molar Tubes feature an intelligent design for more accurate placement, closer fit and stronger bond. The ultra low profile minimizes the potential for debonding due to occlusal forces while improving patient comfort. And the bondable pad with an occlusal indent conforms to the buccal groove for a better fit.



Features

- Extremely low profile for greater patient comfort
- Minimizes potential for occlusal force interference
- Funnel-shaped entrance for ease of wire insertion
- Indent guides for easy, accurate positioning available on select sizes

Find out more at tportho.com
or call **800-348-8856**

*MBT is a registered trademark of 3M Unitek. TP Orthodontics version of the McLaughlin, Bennett and Trevisi prescription is not intended to be a duplication of, nor an endorsement by, these doctors and such endorsement is not implied by the use of this term.

Readi-Base is a registered trademark of TP Orthodontics, Inc. and manufactured under US Patents 6,685,468, 7,131,836, 7,780,442, RE45,743; Britain 1562511; China 03822189.6, ZL 03824889.1; Europe 1545380; France 1562511; Germany 60332263.8-08; Italy 1562511; Japan 4414334, 4503440; Mexico 265725, 274623; Spain 1562511, ES2362089T3. Nu-Edge and eXact are registered trademarks of TP Orthodontics, Inc. NX is a trademark of TP Orthodontics, Inc.

©2021 TP Orthodontics, Inc. All rights reserved.

996-478ENGU, Rev. B



TP Orthodontics, Inc.
tportho.com