

Pre-Finisher®

Semi-custom Finishing Appliance

Efficient finishing for aligner cases

The Pre-Finisher appliance is a semi-custom solution perfect for improving the occlusion in aligner and fixed-appliance patients. When used as an adjunct to aligner treatment, they can reduce or eliminate additional trays required for case refinement – while correcting the posterior open bite often seen with aligner patients.



Nothing could be easier! Near the end of treatment, take a simple measurement and select a Pre-Finisher appliance from your in-office inventory kit – no impressions required. You'll appreciate the time-savings and the final results.

- No laboratory wait or fees; immediate delivery to the patient
- Relates teeth in both arches simultaneously; corrects open bite
- Aligns, rotates and closes spaces
- Sized in one-millimeter increments for precise fitting
- Optional seating springs facilitate placement and retention
- Ideal solution for the already-compliant aligner patient
- Results in 2-3 weeks

Real results



The Refine® Orthodontic Finishing Appliance is also available as a customized solution specially developed for aligner treatment.

Find out more at **tportho.com**
or call **800-348-8856**



TP Orthodontics, Inc.
tportho.com

Pre-Finisher® Semi-custom Finishing Appliance

Parts List

Multi-Appliance Kits

Introductory Kit

(10 Assorted Sizes)	Plain
Nonextraction	190-010

Contents

190-010 Contents	Qty
190-048	1
190-049	1
190-050	2
190-051	1
190-052	2
190-053	1
190-054	1
190-055	1

Inventory Kit

(20 Assorted Sizes)	Plain
Nonextraction	190-020

Contents

190-020 Contents	Qty
190-046	1
190-047	1
190-048	1
190-049	2
190-050	4
190-051	2
190-052	4
190-053	2
190-054	1
190-055	1
190-056	1



Size, in millimeters (last two digits of part number) molded into each appliance.

Individual Appliances

The last two digits of the order number indicate size of Pre-Finisher appliance in millimeters.

Nonextraction Series

Can also be used on first molar extraction cases.

Plain					With Seating Springs				
190-044	190-047	190-050	190-053	190-056	190-144	190-147	190-150	190-153	190-156
190-045	190-048	190-051	190-054	190-057	190-145	190-148	190-151	190-154	190-157
190-046	190-049	190-052	190-055	190-058	190-146	190-149	190-152	190-155	190-158

Fifteen sizes (44mm to 58mm).

Four Bicuspid Extraction Series

Plain				With Seating Springs			
190-449	190-451	190-453	190-455	190-549	190-551	190-553	190-555
190-450	190-452	190-454	190-456	190-550	190-552	190-554	190-556

Eight sizes (49mm to 56mm).

Asian Appliances

Asian appliances are specially designed to address common Asian dental anatomy traits:

- Increased incisal mamelon ridge definition
- Decreased clinical crown length

Asian Nonextraction Series

Plain			
190-848	190-850	190-852	190-854
190-849	190-851	190-853	190-855

Eight sizes (48mm to 55mm).

Asian Four Bicuspid Extraction Series

Plain				
190-644	190-647	190-650	190-653	190-656
190-645	190-648	190-651	190-654	190-657
190-646	190-649	190-652	190-655	

Fourteen sizes (44mm to 57mm).