

TP Orthodontics Indirect Bonding

Precise Bracket Placement

Digitally-designed trays allow for optimal bracket placement.

TPO® Indirect Bonding

Precise Placement. Greater Efficiency.



- Indirect bonding can improve accuracy of bracket placement¹ and reduce the need for repositioning
- End-to-end digital workflow for convenience; 3D-printed trays for accuracy
- Streamlined bonding appointments and reduced chair time
- Brackets are securely positioned, and trays are easy to remove after bonding
- Available on select metal and ceramic brackets in popular prescriptions:
 - InVu and ClearVu aesthetic brackets
 - Nu-Edge twin and self-ligating brackets
 - Tip-Edge metal and ceramic brackets
 - Bondable molar tubes



Try indirect bonding today and see how you can gain a greater level of efficiency in your practice.

Find out more at tportho.com/CUSTOM
or call **800-348-8856**

¹ De Oliveira NS, Rossouw E, Lages EMB, Macari S, Pretti H. "Influence of clinical experience on accuracy of virtual orthodontic attachment bonding in comparison with the direct procedure." Angle Orthod. pp. 1-8: 2019.



TP Orthodontics Indirect Bonding

Precise Bracket Placement

The accuracy you expect,
with the efficiency we deliver.

TPO Indirect Bonding Digital Workflow

- 1 Case Submission**
Submit scan directly from your scanner
Send Rx form to digitalrx@tportho.com
See tportho.com/CUSTOM for details
- 2 TPO Initial Setup**
Our technicians segment the teeth and digitally place brackets
- 3 Doctor Approval via TPRX**
You'll receive notification when your case is ready to review
Approve as-is or request any desired adjustments
- 4 3D Printing**
The tray is printed and brackets are inserted; trays
are sectioned per the Rx form for ease of placement
- 5 Delivery**
Ready to place with your preferred bonding adhesive*



TPRX™ digital custom
appliance management

Submit, manage, and archive
all your prescriptions today.
tportho.com/CUSTOM

Find out more at tportho.com/CUSTOM
or call **800-348-8856**

*Get a free kit of eXact Orthodontic Adhesive to create a most reliable bond.
Available with your first Indirect Bonding case (mention code **EXACT** with
your Rx submission).



TP Orthodontics, Inc.
tportho.com