



TP Orthodontics, Inc.

tportho.com

cosmeticbraces.com

AESTHETIC BRACKET SYSTEMS

AESTHETIC BRACKET SYSTEMS

ClearVu®

- Mini, MBT
- Mini, Roth

ClearVu® with Pre-Coated Adhesive

- Mini, MBT
- Mini, Roth

InVu®

- Mini, MBT
- MBT
- Mini, Roth
- Roth
- Standard Edgewise

InVu® with Pre-Coated Adhesive

- MBT
- Roth

Tip-Edge® PLUS Ceramic

Tip-Edge® PLUS Ceramic with Pre-Coated Adhesive

ClearVu® Mini Ceramic Brackets

MBT MCLAUGHLIN, BENNETT, TREVISI SYSTEM

PATIENT KITS

	With Hooks	0.46mm (.018") Slot	0.56mm (.022") Slot
Kits contain all ceramic brackets			
Upper 3-3	Cuspid	297-380CV	297-381CV
Upper 5-5	Cuspid/Bicuspid	297-408CV	297-409CV
Upper & Lower 3-3	Cuspid	297-382CV	297-383CV
Upper & Lower 5-5	Cuspid/Bicuspid	297-410CV	297-411CV



297-495CV Bi-Dimensional Brackets with Hooks
 0.45 mm (.018") centrals and laterals with 0.56mm (.022") cuspids and bicuspid
 Kit contains 20 brackets.

MAXILLARY BRACKETS

	Slot Size	Left/Right	Bondable
Central	0.46mm (.018")	L	297-310CV
Torque +17°		R	297-311CV
Tip +4°	0.56mm (.022")	L	297-312CV
		R	297-313CV
Lateral	0.46mm (.018")	L	297-314CV
Torque +10°		R	297-315CV
Tip +8°	0.56mm (.022")	L	297-316CV
		R	297-317CV
Cuspid with Distal Hook	0.46mm (.018")	L	297-418CV
Torque 0°		R	297-419CV
Tip +8°	0.56mm (.022")	L	297-420CV
		R	297-421CV
Universal Bicuspid with Hook	0.46mm (.018")	L	297-285CV
Torque -7°		R	297-284CV
Tip 0°	0.56mm (.022")	L	297-287CV
In/Out Comp. .027"		R	297-286CV

Package of 5

MANDIBULAR BRACKETS

	Slot Size	Left/Right	Bondable
Anterior Narrow	0.46mm (.018")	L/R	297-354CV
Torque -6°			
Tip 0°	0.56mm (.022")	L/R	297-355CV
Cuspid with Distal Hook	0.46mm (.018")	L	297-446CV
Torque 0°		R	297-447CV
Tip +3°	0.56mm (.022")	L	297-448CV
		R	297-449CV
1st Bicuspid with Hook	0.46mm (.018")	L	297-296CV
Torque -12°		R	297-297CV
Tip 0°	0.56mm (.022")	L	297-298CV
		R	297-299CV
2nd Bicuspid with Hook	0.46mm (.018")	L	297-288CV
Torque -17°		R	297-289CV
Tip 0°	0.56mm (.022")	L	297-290CV
		R	297-291CV

Package of 5

PRE-COATED PATIENT KITS

	With Hooks	0.46mm (.018") Slot	0.56mm (.022") Slot
Kits contain all ceramic brackets			
Upper & Lower 5-5	Cuspid/Bicuspid	297-410CVX	297-411CVX



MBT is a registered trademark of 3M. TP Orthodontics version of the McLaughlin, Bennett and Trevisi prescription is not intended to be a duplication of, nor an endorsement by, these doctors and such endorsement is not implied by the use of this term.

800-348-8856

ClearVu® Mini Ceramic Brackets

ROTH SYSTEM

PATIENT KITS

	With Hooks	0.46mm (.018") Slot	0.56mm (.022") Slot
Kits contain all ceramic brackets			
Upper 3-3	Cuspid	297-100CV	297-101CV
Upper 5-5	Cuspid/Bicuspid	297-404CV	297-405CV
Upper & Lower 3-3	Cuspid	297-102CV	297-103CV
Upper & Lower 5-5	Cuspid/Bicuspid	297-406CV	297-407CV



297-490CV Bi-Dimensional Brackets with Hooks
 0.46mm (.018") centrals and laterals with 0.56mm (.022") cuspids and bicuspid
 Kit contains 20 brackets.

MAXILLARY BRACKETS

	Slot Size	Left/Right	Bondable
Central	0.46mm (.018")	L	297-138CV
Torque +14°		R	297-139CV
Tip +5°	0.56mm (.022")	L	297-150CV
		R	297-151CV
Lateral	0.46mm (.018")	L	297-148CV
Torque +7°		R	297-149CV
Tip +8°	0.56mm (.022")	L	297-152CV
		R	297-153CV
Cuspid with Distal Hook	0.46mm (.018")	L	297-418CV
Torque 0°		R	297-419CV
Tip +8°	0.56mm (.022")	L	297-420CV
		R	297-421CV
Universal Bicuspid with Hook	0.46mm (.018")	L	297-284CV
Torque -7°		R	297-285CV
Tip 0°	0.56mm (.022")	L	297-286CV
		R	297-287CV

Package of 5

MANDIBULAR BRACKETS

	Slot Size	Left/Right	Bondable
Anterior	0.46mm (.018")	L/R	297-169CV
Torque -1°			
Tip 0°	0.56mm (.022")	L/R	297-160CV
Cuspid with Distal Hook	0.46mm (.018")	L	297-426CV
Torque -11°		R	297-427CV
Tip +6°	0.56mm (.022")	L	297-428CV
		R	297-429CV
1st Bicuspid with Hook	0.46mm (.018")	L	297-288CV
Torque -17°		R	297-289CV
Tip 0°	0.56mm (.022")	L	297-290CV
		R	297-291CV
2nd Bicuspid with Hook	0.46mm (.018")	L	297-292CV
Torque -22°		R	297-293CV
Tip 0°	0.56mm (.022")	L	297-294CV
		R	297-295CV

Package of 5

PRE-COATED PATIENT KITS

	With Hooks	0.46mm (.018") Slot	0.56mm (.022") Slot
Kits contain all ceramic brackets			
Upper & Lower 5-5	Cuspid/Bicuspid	297-406CVX	297-407CVX
Upper & Lower 5-5	Cuspid		297-414CVX




Compatible with the Roth system: not endorsed by Dr. Roth.


ClearVu® Ceramic Brackets

Bracket Identification Charts

MBT



Roth



MBT is a registered trademark of 3M. TP Orthodontics version of the McLaughlin, Bennett and Trevisi prescription is not intended to be a duplication of, nor an endorsement by, these doctors and such endorsement is not implied by the use of this term.

Compatible with the Roth system: not endorsed by Dr. Roth.

800-348-8856

InVu® Mini Ceramic Brackets

MBT MCLAUGHLIN, BENNETT, TREVISI SYSTEM

PATIENT KITS

	With Hooks	0.46mm (.018") Slot	0.56mm (.022") Slot
Kits contain all ceramic brackets			
Upper 3-3	No Hooks	297-391	297-393
Upper 5-5	Cuspid	297-394	297-396
Upper 5-5	Cuspid/Bicuspid	297-408	297-409
Upper & Lower 3-3	No Hooks	297-303	297-304
Upper & Lower 5-5	Cuspid	297-306	297-308
Upper & Lower 5-5	Cuspid/Bicuspid	297-410	297-411
Upper & Lower 5-5	No Hooks	297-307	297-309



MAXILLARY BRACKETS

	Slot Size	Left/Right	Bondable
Central	0.46mm	L	297-310
Torque +17°	(.018")	R	297-311
Tip +4°	0.56mm	L	297-312
	(.022")	R	297-313
Lateral	0.46mm	L	297-314
Torque +10°	(.018")	R	297-315
Tip +8°	0.56mm	L	297-316
	(.022")	R	297-317
Cuspid without Hook	0.46mm	L	297-318
Torque -7°	(.018")	R	297-319
Tip +8°	0.56mm	L	297-320
	(.022")	R	297-321
Cuspid with Distal Hook	0.46mm	L	297-418
Torque 0°	(.018")	R	297-419
Tip +8°	0.56mm	L	297-420
	(.022")	R	297-421
Universal Bicuspid w/o Hook	0.46mm	L/R	297-226M
Torque -7°	(.018")		
Tip 0°	0.56mm	L/R	297-228M
	(.022")		
Universal Bicuspid w Hook	0.46mm	L	297-285M
Torque -7°	(.018")	R	297-284M
Tip 0°	0.56mm	L	297-287M
	(.022")	R	297-286M

Package of 5

MANDIBULAR BRACKETS

	Slot Size	Left/Right	Bondable
Anterior	0.46mm	L/R	297-354
Torque -6°	(.018")		
Tip 0°	0.56mm	L/R	297-355
	(.022")		
Cuspid without Hook	0.46mm	L	297-356
Torque -6°	(.018")	R	297-357
Tip +3°	0.56mm	L	297-358
	(.022")	R	297-359
Cuspid with Distal Hook	0.46mm	L	297-446
Torque 0°	(.018")	R	297-447
Tip +3°	0.56mm	L	297-448
	(.022")	R	297-449
1st Bicuspid without Hook	0.46mm	L/R	297-264M
Torque -12°	(.018")		
Tip 0°	0.56mm	L/R	297-266M
	(.022")		
1st Bicuspid with Hook	0.46mm	L	297-296M
Torque -12°	(.018")	R	297-297M
Tip 0°	0.56mm	L	297-298M
	(.022")	R	297-299M
2nd Bicuspid without Hook	0.46mm	L/R	297-272M
Torque -17°	(.018")		
Tip 0°	0.56mm	L/R	297-274M
	(.022")		
2nd Bicuspid with Hook	0.46mm	L	297-288M
Torque -17°	(.018")	R	297-289M
Tip 0°	0.56mm	L	297-290M
	(.022")	R	297-291M


Package of 5

InVu® Ceramic Brackets

MBT MCLAUGHLIN, BENNETT, TREVISI SYSTEM

PATIENT KITS

	With Hooks	0.46mm (.018") Slot	0.56mm (.022") Slot
Kits contain all ceramic brackets			
Upper 3-3	Cuspid	297-190	297-192
	No Hooks		297-193
Upper 5-5	Cuspid	297-194	297-196
	Cuspid/Bicuspid	297-400	297-401
	No Hooks		297-197
Upper & Lower 3-3	Cuspid	297-202	297-204
Upper & Lower 5-5	Cuspid	297-206	297-208
	Cuspid/Bicuspid	297-402	297-403



MAXILLARY BRACKETS

	Slot Size	Left/Right	Bondable
Central	0.46mm (.018")	L	297-210
		R	297-211
	0.56mm (.022")	L	297-212
Lateral	0.46mm (.018")	L	297-214
		R	297-215
	0.56mm (.022")	L	297-216
Cuspid without Hook	0.46mm (.018")	L	297-218
		R	297-219
	0.56mm (.022")	L	297-220
Cuspid with Distal Hook	0.46mm (.018")	L	297-222
		R	297-223
	0.56mm (.022")	L	297-224
Universal Bicuspid without Hook	0.46mm (.018")	L/R	297-226
		L/R	297-228
	0.56mm (.022")	L/R	297-228
Universal Bicuspid with Hook	0.46mm (.018")	L	297-285
		R	297-284
	0.56mm (.022")	L	297-287
In/Out Comp. .027"	0.56mm (.022")	R	297-286

Package of 5

MANDIBULAR BRACKETS

	Slot Size	Left/Right	Bondable
Anterior Narrow	0.46mm (.018")	L/R	297-254
		L/R	297-255
	0.56mm (.022")	L/R	297-255
Cuspid without Hook	0.46mm (.018")	L	297-256
		R	297-257
	0.56mm (.022")	L	297-258
Cuspid with Distal Hook	0.46mm (.018")	L	297-260
		R	297-261
	0.56mm (.022")	L	297-262
1st Bicuspid without Hook	0.46mm (.018")	L/R	297-264
		L/R	297-266
	0.56mm (.022")	L/R	297-266
1st Bicuspid with Hook	0.46mm (.018")	L	297-296
		R	297-297
	0.56mm (.022")	L	297-298
2nd Bicuspid without Hook	0.46mm (.018")	L/R	297-272
		L/R	297-274
	0.56mm (.022")	L/R	297-274
2nd Bicuspid with Hook	0.46mm (.018")	L	297-288
		R	297-289
	0.56mm (.022")	L	297-290
In/Out Comp. .027"	0.56mm (.022")	R	297-291

Package of 5

MBT is a registered trademark of 3M. TP Orthodontics version of the McLaughlin, Bennett and Trevisi prescription is not intended to be a duplication of, nor an endorsement by, these doctors and such endorsement is not implied by the use of this term.


800-348-8856

InVu® Mini Ceramic Brackets

ROTH SYSTEM

PATIENT KITS

	With Hooks	0.46mm (.018") Slot	0.56mm (.022") Slot
Kits contain all ceramic brackets			
Upper 3-3	No Hooks	297-106	297-110
Upper 4-4	Cuspid		297-057
Upper 5-5	Cuspid	297-060	297-065
Upper 5-5	Cuspid/Bicuspid	297-404	297-405
Upper 5-5	No Hooks	297-107	297-111
Upper & Lower 3-3	No Hooks	297-108	297-112
Upper & Lower 5-5	Cuspid	297-070	297-075
Upper & Lower 5-5	Cuspid/Bicuspid	297-406	297-407
Upper & Lower 5-5	No Hooks	297-109	297-113



MAXILLARY BRACKETS

	Slot Size	Left/Right	Bondable
Central	0.46mm	L	297-138
Torque +14°	(.018")	R	297-139
Tip +5°	0.56mm	L	297-150
	(.022")	R	297-151
Lateral	0.46mm	L	297-148
Torque +7°	(.018")	R	297-149
Tip +8°	0.56mm	L	297-152
	(.022")	R	297-153
Cuspid without Hook	0.46mm	L	297-158
Torque 0°	(.018")	R	297-159
Tip +8°	0.56mm	L	297-154
	(.022")	R	297-155
Cuspid with Distal Hook	0.46mm	L	297-418
Torque 0°	(.018")	R	297-419
Tip +8°	0.56mm	L	297-420
	(.022")	R	297-421
Universal Bicuspid w/o Hook	0.46mm	L/R	297-226M
Torque -7°	(.018")		
Tip 0°	0.56mm	L/R	297-228M
	(.022")		
Universal Bicuspid w/ Hook	0.46mm	L	297-284M
Torque -7°	(.018")	R	297-285M
Tip 0°	0.56mm	L	297-286M
	(.022")	R	297-287M

Package of 5

MANDIBULAR BRACKETS

	Slot Size	Left/Right	Bondable
Anterior	0.46mm	L/R	297-169
Torque -1°	(.018")		
Tip 0°	0.56mm	L/R	297-160
	(.022")		
Cuspid without Hook	0.46mm	L	297-176
Torque -11°	(.018")	R	297-177
Tip +6°	0.56mm	L	297-174
	(.022")	R	297-175
Cuspid with Distal Hook	0.46mm	L	297-426
Torque -11°	(.018")	R	297-427
Tip +6°	0.56mm	L	297-428
	(.022")	R	297-429
1st Bicuspid without Hook	0.46mm	L/R	297-182M
Torque -17°	(.018")		
Tip 0°	0.56mm	L/R	297-184M
	(.022")		
1st Bicuspid with Hook	0.46mm	L	297-288M
Torque -17°	(.018")	R	297-289M
Tip 0°	0.56mm	L	297-290M
	(.022")	R	297-291M
2nd Bicuspid without Hook	0.46mm	L/R	297-186M
Torque -22°	(.018")		
Tip 0°	0.56mm	L/R	297-188M
	(.022")		
2nd Bicuspid with Hook	0.46mm	L	297-292M
Torque -22°	(.018")	R	297-293M
Tip 0°	0.56mm	L	297-294M
	(.022")	R	297-295M


Package of 5

InVu® Ceramic Brackets

ROTH SYSTEM

PATIENT KITS

	With Hooks	0.46mm (.018") Slot	0.56mm (.022") Slot
Kits contain all ceramic brackets			
Upper 3-3	Cuspid	297-120	297-121
Upper 4-4	Cuspid	297-122	
Upper 5-5	Cuspid	297-124	297-125
Upper 5-5	Cuspid/Bicuspid	297-126	297-127
Upper & Lower 3-3	Cuspid	297-132	297-133
Upper & Lower 5-5	Cuspid	297-136	297-137
Upper & Lower 5-5	Cuspid/Bicuspid	297-128	297-129



MAXILLARY BRACKETS

	Slot Size	Left/Right	Bondable
Central	0.46mm	L	297-140
Torque +14°	(.018")	R	297-141
Tip +5°	0.56mm	L	297-142
	(.022")	R	297-143
Lateral	0.46mm	L	297-144
Torque +7°	(.018")	R	297-145
Tip +8°	0.56mm	L	297-146
	(.022")	R	297-147
Cuspid with Distal Hook	0.46mm	L	297-222
Torque 0°	(.018")	R	297-223
Tip +8°	0.56mm	L	297-224
	(.022")	R	297-225
Universal Bicuspid w/o Hook	0.46mm	L/R	297-226
Torque -7°	(.018")		
Tip 0°	0.56mm	L/R	297-228
	(.022")		
Universal Bicuspid with Hook	0.46mm	L	297-284
Torque -7°	(.018")	R	297-285
Tip 0°	0.56mm	L	297-286
	(.022")	R	297-287

Package of 5

MANDIBULAR BRACKETS

	Slot Size	Left/Right	Bondable
Anterior	0.46mm	L/R	297-170
Torque -1°	(.018")		
Tip 0°	0.56mm	L/R	297-172
	(.022")		
Cuspid with Distal Hook	0.46mm	L	297-178
Torque -11°	(.018")	R	297-179
Tip +6°	0.56mm	L	297-180
	(.022")	R	297-181
1st Bicuspid without Hook	0.46mm	L/R	297-182
Torque -17°	(.018")		
Tip 0°	0.56mm	L/R	297-184
	(.022")		
1st Bicuspid with Hook	0.46mm	L	297-288
Torque -17°	(.018")	R	297-289
Tip 0°	0.56mm	L	297-290
	(.022")	R	297-291
2nd Bicuspid without Hook	0.46mm	L/R	297-186
Torque -22°	(.018")		
Tip 0°	0.56mm	L/R	297-188
	(.022")		
2nd Bicuspid with Hook	0.46mm	L	297-292
Torque -22°	(.018")	R	297-293
Tip 0°	0.56mm	L	297-294
	(.022")	R	297-295

Package of 5

Compatible with the Roth system: not endorsed by Dr. Roth.


800-348-8856

InVu® Ceramic Brackets

STANDARD EDGEWISE SYSTEM

PATIENT KITS

	0.46mm (.018") Slot	0.56mm (.022") Slot
Kits contain all ceramic brackets without cuspid hooks		
Upper 3-3	297-960	297-961
Upper & Lower 5-5	297-964	297-965



MAXILLARY BRACKETS

		Slot Size	Left/ Right	Bondable
Central		0.46mm (.018")	L/R	297-905
Torque	0°			
Tip	0°	0.56mm (.022")	L/R	297-906
Lateral		0.46mm (.018")	L/R	297-911
Torque	0°			
Tip	0°	0.56mm (.022")	L/R	297-912
Cuspid		0.46mm (.018")	L/R	297-929
Torque	0°			
Tip	0°	0.56mm (.022")	L/R	297-930
Bicuspid without Hook		0.46mm (.018")	L/R	297-939
Torque	0°			
Tip	0°	0.56mm (.022")	L/R	297-940

Package of 5

MANDIBULAR BRACKETS

		Slot Size	Left/ Right	Bondable
Anterior		0.46mm (.018")	L/R	297-933
Torque	0°			
Tip	0°	0.56mm (.022")	L/R	297-934
Cuspid		0.46mm (.018")	L/R	297-929
Torque	0°			
Tip	0°	0.56mm (.022")	L/R	297-930
Bicuspid without Hook		0.46mm (.018")	L/R	297-941
Torque	0°			
Tip	0°	0.56mm (.022")	L/R	297-942


Package of 5

InVu® Ceramic Brackets

MBT MCLAUGHLIN, BENNETT, TREVISI SYSTEM

PRE-COATED PATIENT KITS — MINI

	With Hooks	0.46mm (.018") Slot	0.56mm (.022") Slot
Kits contain all ceramic brackets			
Upper 5-5	Cuspid/Bicuspid	297-408RBX	
Upper & Lower 5-5	Cuspid		297-308RBX



PRE-COATED PATIENT KITS

	With Hooks	0.46mm (.018") Slot	0.56mm (.022") Slot
Kits contain all ceramic brackets			
Upper 5-5	Cuspid/Bicuspid	297-400RBX	297-401RBX
Upper & Lower 5-5	Cuspid/Bicuspid	297-402RBX	297-403RBX

InVu® Ceramic Brackets

ROTH SYSTEM

PRE-COATED PATIENT KITS - MINI

	With Hooks	0.46mm (.018") Slot
Kits contain all ceramic brackets		
Upper 5-5	Cuspid/Bicuspid	297-404RBX
Upper & Lower 5-5	Cuspid	297-070RBX



PRE-COATED PATIENT KITS

	With Hooks	0.46mm (.018") Slot	0.56mm (.022") Slot
Kits contain all ceramic brackets			
Upper 5-5	Cuspid/Bicuspid	297-126RBX	
Upper & Lower 5-5	Cuspid	297-136RBX	297-137RBX

MBT is a registered trademark of 3M. TP Orthodontics version of the McLaughlin, Bennett and Trevisi prescription is not intended to be a duplication of, nor an endorsement by, these doctors and such endorsement is not implied by the use of this term.


Compatible with the Roth system: not endorsed by Dr. Roth.

800-348-8856


InVu® Ceramic Brackets

Bracket Identification Charts


MBT



Roth



Standard Edgewise



MBT is a registered trademark of 3M. TP Orthodontics version of the McLaughlin, Bennett and Trevisi prescription is not intended to be a duplication of, nor an endorsement by, these doctors and such endorsement is not implied by the use of this term.

Compatible with the Roth system: not endorsed by Dr. Roth.

Tip-Edge® PLUS Ceramic Brackets



PATIENT KITS

	Slot Size	Non-Oriented	Oriented Incisal Jig
Kits contain all ceramic brackets			
Upper & Lower 5-5	0.56mm		297-602
Upper 3-3	(.022")		297-598
Upper & Lower 3-3			297-599
Upper 5-5 Non-extraction			297-600
Lower 3-3			297-604
Lower 5-5 Non-extraction			297-605

	Slot Size	Non-Oriented	Oriented Incisal Jig
Kits combine ceramic with metal brackets			
Up. 3-3 Ceramic & Up. 4-5 Metal with Lo. 5-5 Metal			
Non-extraction	0.56mm	297-608	297-610
	(.022")		
Up. 5-5 Ceramic with Lo. 5-5 Metal			
Non-extraction	0.56mm		297-614
	(.022")		

Non-extraction kits contain 20 bondable brackets, extraction kits have 16.

MAXILLARY BRACKETS

	Slot Size	Left/Right	Oriented 3.5mm Blue	Oriented 4mm Red	Oriented 4.5mm Yellow	Oriented 5mm Green
Central						
Torque +12°	0.56mm	L		297-639	297-655	
Tip +5°	(.022")	R		297-640	297-656	
Lateral						
Torque +8°	0.56mm	L	297-629	297-641		
Tip +9°	(.022")	R	297-630	297-642		
Cuspid						
Torque -4°	0.56mm	L		297-643	297-659	297-671
Tip +11°	(.022")	R		297-644	297-660	297-672
Bicuspid						
Torque -7°	0.56mm	CCW	297-631	297-645	297-661	
Tip 0°	(.022")	CW	297-632	297-646	297-662	

Package of 5

MANDIBULAR BRACKETS

Anterior						
Torque -1°	0.56mm	L	297-633	297-647		
Tip +2°	(.022")	R	297-634	297-648		
Cuspid						
Torque -11°	0.56mm	L		297-651	297-667	297-673
Tip +5°	(.022")	R		297-652	297-668	297-674
Bicuspid						
Torque -20°	0.56mm	CCW	297-637	297-653	297-669	
Tip 0°	(.022")	CW	297-638	297-654	297-670	

Recommended jig heights: Maxillary central 4mm, lateral 3.5mm, cuspid 4.5mm, bicuspid 4mm.

Mandibular anteriors 3.5mm, cuspids 4mm, bicuspid 4mm.

Package of 5

PRE-COATED PATIENT KIT


	Slot Size	Oriented Incisal Jig
Kits contain all ceramic brackets		
Upper & Lower 5-5	0.56mm	RB297-602
Non-extraction	(.022")	

Pre-coated
with eXact®
Adhesive


Tip-Edge® PLUS Ceramic Brackets

Bracket Identification Charts

Non-Extraction



Extraction




Circles identify the upper brackets and triangles identify the lower brackets. Arrows on the bicuspids represent the direction crowns should tip during initial treatment – either clockwise or counterclockwise. The bonding pad shape also points toward the direction of desired crown tipping.


Positioning Jigs

Incisal jigs are designed to help with accurate incisogingival bracket placement and assist with sighting long axis crown angulations.

Jig measurements are from the center of the archwire slot to the incisal edge or cusp tips. Color coding identifies the various heights.



RECOMMENDED JIG HEIGHT



Jig height substitutions can be made for an additional charge. Green 5 mm jigs are also available for cuspids.

Gingivally offset mandibular bicuspid brackets are available with brackets positioned gingivally on bonding pads to permit proper bonding height on partially erupted teeth