



TP Orthodontics, Inc.

tportho.com

cosmeticbraces.com

TUBES, ATTACHMENTS AND BANDS

TUBES, ATTACHMENTS AND BANDS

Ultra Low Profile Single Molar Tubes

- First Molar
- Second Molar

Standard Single Molar Tubes

- First Molar
- Convertible
- Second Molar

Headgear and Auxiliary Molar Tubes

- Combination Convertible
- Combination Convertible Rotation
- Double Convertible
- Combination Convertible Lip Bumper
- Triple Convertible Rotation
- Triple .036" Tip-Edge Auxiliary

Tip-Edge® System Molar Tubes

- Combination
- Combination Convertible

Lingual Attachments

- Bondable
- Weldable

GripTite® First Molar Bands

Pre-Fit® Bands

- Second Molar
- Bicuspid

Prewelding Service

Prewelded Band Assemblies

Ultra Low Profile Single Molar Tubes

First Molar

6 | 6 UPPER FIRST MOLAR—ULTRA LOW PROFILE

Archwire Tube Size	Tube Length	Torque	Tip	Distal Offset	Left/ Right	Weldable*	Bondable
0.46 x 0.64mm (.018" x .025")	4.0	-10°	0°	0°	L	391-738W	391-742
					R	391-739W	391-743
0.56 x 0.71mm (.022" x .028")	4.0	-10°	0°	0°	L	391-740W	391-746
					R	391-741W	391-747
	4.5	-14°	0°	10°	L		391-936
					R		391-937

Package of 10

*Weldable (W) tubes are only available individually welded to TPO bands.



First Molar Tube

6 | 6 LOWER FIRST MOLAR—ULTRA LOW PROFILE

Archwire Tube Size	Tube Length	Torque	Tip	Distal Offset	Left/ Right	Weldable*	Bondable
0.46 x 0.64mm (.018" x .025")	4.0	-10°	0°	0°	L	391-739W	391-748
					R	391-738W	391-749
		-20°	0°	0°	L		391-816
					R		391-817
0.56 x 0.71mm (.022" x .028")	4.0	-10°	0°	0°	L	391-741W	391-750
					R	391-740W	391-751
		-20°	0°	0°	L	391-820W	391-822
					R	391-821W	391-823
	4.5	-25°	0°	0°	L		391-944
					R		391-945

Package of 10

*Weldable (W) tubes are only available individually welded to TPO bands.

Second Molar

7 | 7 UPPER SECOND MOLAR—ULTRA LOW PROFILE

Archwire Tube Size	Tube Length	Torque	Tip	Distal Offset	Left/ Right	Bondable
0.56 x 0.71mm (.022" x .028")	3.2	-14°	0°	10°	L	391-952
					R	391-953

Package of 10



Second Molar Tube

7 | 7 LOWER SECOND MOLAR—ULTRA LOW PROFILE

Archwire Tube Size	Tube Length	Torque	Tip	Distal Offset	Left/ Right	Bondable
0.56 x 0.71mm (.022" x .028")	3.2	-25°	0°	5°	L	391-960
					R	391-961

Package of 10

Standard Single Molar Tubes

First Molar

6|6 UPPER FIRST MOLAR WITH DISTAL EXTENSION

Archwire Tube Size	Tube Length	Torque	Tip	Distal Offset	Left/ Right	Weldable*	Bondable	
							Original Pad .182"	Big Foot Pad .325"
0.46 x 0.64mm (.018" x .025")	4.2	-10°	0°	14°	L		391-600	391-490
		-14°	0°	14°	R		391-601	391-491
0.56 x 0.71mm (.022" x .028")	4.2	-10°	0°	14°	L	391-624W	391-616	391-494
		-14°	0°	14°	R	391-625W	391-617	391-495
					L			391-496
					R			391-497

Package of 10

*Weldable (W) tubes are only available individually welded to TPO bands.

Upper Right



6|6 LOWER FIRST MOLAR WITH DISTAL EXTENSION

Archwire Tube Size	Tube Length	Torque	Tip	Distal Offset	Left/ Right	Weldable*	Bondable	
							Original Pad .182"	Big Foot Pad .325"
0.46 x 0.64mm (.018" x .025")	4.2	-25°	0°	4°	L			391-498
		-30°	2°	4°	R	391-614W		391-499
0.56 x 0.71mm (.022" x .028")	4.2	-25°	0°	4°	L	391-628W	391-620	391-578
					R	391-629W	391-621	391-579

Package of 10

*Weldable (W) tubes are only available individually welded to TPO bands.

Lower Right



6|6 UPPER FIRST MOLAR

Archwire Tube Size	Tube Length	Torque	Tip	Distal Offset	Left/ Right	Weldable*	Bondable
							Original Pad .260"
0.46 x 0.64mm (.018" x .025")	4.0	0°	0°	0°	L		210-396L
					R		210-396R
0.56 x 0.71mm (.022" x .028")	4.0	0°	0°	0°	L	391-012WL	210-397L
					R	391-012WR	210-397R

Package of 10

*Weldable (W) tubes are only available individually welded to TPO bands.

Upper Right



6|6 LOWER FIRST MOLAR

Archwire Tube Size	Tube Length	Torque	Distal Tip	Offset	Left/ Right	Weldable*	Bondable
							Original Pad .260"
0.46 x 0.64mm (.018" x .025")	4.0	0°	0°	0°	L		210-396R
					R		210-396L
0.56 x 0.71mm (.022" x .028")	4.0	0°	0°	0°	L	391-012WR	210-397R
					R	391-012WL	210-397L

Package of 10

*Weldable (W) tubes are only available individually welded to TPO bands.

Lower Right



Standard Single Molar Tubes Convertible

6/6 MICRO6 COBALT CHROMIUM UPPER CONVERTIBLE FIRST MOLAR WITH ZIP-RELEASE

Archwire Tube Size	Tube Length	Torque	Tip	Distal Offset	Left/ Right	Bondable
						Standard Pad
0.46 x 0.64mm (.018" x .025")	4.0	-10°	0°	4°	L	391-980
					R	391-981
0.56 x 0.71mm (.022" x .028")	4.0	-10°	0°	4°	L	391-984
					R	391-985

Package of 10

Upper Right



6/6 MICRO6 COBALT CHROMIUM UPPER CONVERTIBLE FIRST MOLAR

Archwire Tube Size	Tube Length	Torque	Distal Tip	Left/ Offset	Right	Bondable		
						Standard Pad	Small Foot Pad .260"	Big Foot Pad .325"
0.46 x 0.64mm (.018" x .025")	4.0	-10°	0°	4°	L	391-528	391-652	391-660
					R	391-529	391-653	391-661
0.56 x 0.71mm (.022" x .028")	4.0	-10°	0°	4°	L	391-701	391-654	391-662
					R	391-702	391-655	391-663

Package of 10

Upper Right



6/6 LOWER CONVERTIBLE FIRST MOLAR

Archwire Tube Size	Tube Length	Torque	Tip	Distal Offset	Left/ Right	Weldable*	Bondable	
							Original Pad .182"	Big Foot Pad .325"
0.46 x 0.64mm (.018" x .025")	4.2	-25°	0°	4°	L	391-536W		391-664
					R	391-537W		391-665
0.56 x 0.71mm (.022" x .028")	4.2	-25°	0°	4°	L	391-540W	391-548	391-666
					R	391-541W	391-549	391-667

Package of 10

*Weldable (W) tubes are only available individually welded to TPO bands.

Lower Right



Standard Single Molar Tubes

Second Molar

7/7 UPPER AND LOWER SECOND MOLARS MINI SIZE

Upper/ Lower	Archwire Tube Size	Tube Length	Torque	Tip	Distal Offset	Left/ Right	Weldable*	Bondable
Upper	0.46 x 0.64mm (.018" x .025")	2.3	-14°	0°	10°	L		391-632
						R		391-633
	0.56 x 0.71mm (.022" x .028")	2.3	-10°	0°	6°	L		391-784
						R		391-785
Lower	0.46 x 0.64mm (.018" x .025")	2.3	-14°	0°	10°	L	391-794W	391-634
						R	391-795W	391-635
	0.56 x 0.71mm (.022" x .028")	2.3	-10°	0°	0°	L		391-636
						R		391-637
			-10°	0°	0°	L	391-798W	391-638
						R	391-799W	391-639
			-25°	0°	6°	L		391-886
						R		391-887

Package of 10

*Weldable (W) tubes are only available individually welded to TPO bands.

Right



7/7 UPPER SECOND MOLAR WITHOUT DISTAL EXTENSION

Archwire Tube Size	Tube Length	Torque	Tip	Distal Offset	Left/ Right	Weldable*	Bondable
0.46 x 0.64mm (.018" x .025")	3.3	-10°	0°	14°	L		391-394
					R		391-393
0.56 x 0.71mm (.022" x .028")	3.3	-10°	0°	14°	L	391-400W	391-402
					R	391-399W	391-401
		-14°	0°	14°	L		391-406
					R		391-405

Package of 10

*Weldable (W) tubes are only available individually welded to TPO bands.

Upper Right



7/7 LOWER SECOND MOLAR WITHOUT DISTAL EXTENSION

Archwire Tube Size	Tube Length	Torque	Tip	Distal Offset	Left/ Right	Weldable*	Bondable
0.46 x 0.64mm (.018" x .025")	3.3	-25°	0°	4°	L		391-176
					R		391-175
0.56 x 0.71mm (.022" x .028")	3.3	-25°	0°	4°	L	391-138W	391-140
					R	391-137W	391-139
		-30°	0°	4°	L	391-218W	391-220
					R	391-217W	391-219

Package of 10

*Weldable (W) tubes are only available individually welded to TPO bands.

Lower Right



Headgear and Auxiliary Molar Tubes Combination Convertible


6/6 UPPER FIRST MOLAR MINI SIZE OCCLUSAL HEADGEAR TUBE

Main Tube Size	Headgear Tube Size	Torque	Tip	Distal Offset		Left/ Right	Weldable*
				Main	Headgear		
.018" x .025" x 3.5mm	.045" x 3.5mm	-10°	0°	7°	0°	L	391-502W
						R	391-503W
.022" x .028" x 3.5mm	.045" x 3.5mm	-10°	0°	14°	0°	L	391-508W
						R	391-509W

*Weldable (W) tubes are only available individually welded to TPO bands.

Recommended for use with VektorPRO Class II Corrector.

Upper Right



6/6 UPPER FIRST MOLAR MINI SIZE GINGIVAL HEADGEAR TUBE

Main Tube Size	Headgear Tube Size	Torque	Tip	Distal Offset		Left/ Right	Weldable*
				Main	Headgear		
.022" x .028" x 3.5mm	.045" x 3.5mm	-10°	0°	14°	0°	L	391-268W
						R	391-269W

Package of 10

*Weldable (W) tubes are only available individually welded to TPO bands.

Upper Right



Combination Convertible Rotation

6/6 UPPER FIRST MOLAR MINI SIZE ROTATES MOLARS WITHOUT ARCH BARS



Occlusal headgear tube has 7° offset at end of treatment.

Main Tube Size	Headgear Tube Size	Torque	Tip	Distal Offset		Left/ Right	Weldable*
				Main	Headgear		
.018" x .025" x 3.5mm	.045" x 3.5mm	-14°	0°	14°	7°	L	391-510
						R	391-511
		-10°	0°	7°	7°	L	391-512
						R	391-513
.022" x .028" x 3.5mm	.045" x 3.5mm	-14°	0°	14°	7°	L	391-514
						R	391-515
		-10°	0°	7°	7°	L	391-516
						R	391-517

Package of 10

*Can be ordered on TPO bands; add W to the part number.

Upper Right



Auxiliary Molar Tubes


Double Convertible

6|6 UPPER FIRST MOLAR

Main Occl. Tube Size	Aux. Ging. Tube Size	Torque		Tip		Distal Offset		Left/ Right	Weldable*
		Main	Aux.	Main	Aux.	Main	Aux.		
.018" x .025" x 3.5mm	.018" x .025" x 3.5mm	-14°	0°	0°	0°	10°	0°	L	391-680W
								R	391-681W
.022" x .028" x 3.5mm	.022" x .028" x 3.5mm	-14°	0°	0°	0°	10°	0°	L	391-682W
								R	391-683W

*Weldable (W) tubes are only available individually welded to TPO bands.

Upper Right



6|6 LOWER FIRST MOLAR WITH DISTAL EXTENSION WITH VERTICAL SLOT

Main Occl. Tube Size	Aux. Ging. Tube Size	Torque		Tip		Distal Offset		Left/ Right	Weldable*
		Main	Aux.	Main	Aux.	Main	Aux.		
.018" x .025" x 4.2mm	.018" x .025" x 4.2mm	-27°	-15°	0°	0°	5°	0°	L	391-206W
								R	391-205W
.022" x .028" x 4.2mm	.022" x .028" x 4.2mm	-27°	-15°	0°	0°	5°	0°	L	391-648W
								R	391-649W

*Weldable (W) tubes are only available individually welded to TPO bands.

Lower Right



Headgear and Auxiliary Molar Tubes

Combination Convertible Lip Bumper

6/6 LOWER FIRST MOLAR AUTOMATIC ROTATION WITHOUT LIP BUMPER REACTIVATION

- Lip bumper tube has elliptical opening for ease of wire insertion on rotated molars.
- Stops on lip bumper must be larger than openings in mesial ends of tubes.

Lower Right



Main Tube Size	Headgear Tube Size	Torque	Tip	Distal Offset		Left/Right	Weldable*
				Main	Headgear		
.018" x .025" x 4.4mm	.045" x 4.4mm	-27°	0°	4°	0°	L	391-270
						R	391-271
.022" x .028" x 4.4mm	.045" x 4.4mm	-27°	0°	4°	0°	L	391-272
						R	391-273

Package of 10


*Can be ordered on TPO bands; add W to the part number.

Triple Convertible Rotation

6/6 UPPER FIRST MOLAR ROTATES MOLAR WITHOUT ARCH BARS

- Unique elliptical mesial opening of headgear tube facilitates engagement, especially when molars are rotated.
- Eliminates use of palatal arch bars and time-consuming adjustments.
- Headgear tube has 7° distal offset at end of treatment.


Upper Right



Main Tube Size	Aux. Ging. Tube Size	Headgear Tube Size	Torque		Tip		Distal Offset		Left/Right	Weldable*
			Main	Aux.	Main	Aux.	Main	Aux.		
.018" x .025" x 3.5mm	.018" x .025" x 3.5mm	.045" x 3.5mm	-14°	0°	0°	0°	14°	0°	L	391-568
									R	391-569
.022" x .028" x 3.5mm	.022" x .028" x 3.5mm	.045" x 3.5mm	-14°	0°	0°	0°	14°	0°	L	391-432
									R	391-433

Package of 10

*Can be ordered on TPO bands; add W to the part number.



Triple .036" Tip-Edge Auxiliary with Easy-Out Feature

6/6 UPPER FIRST MOLAR

Upper Right



Main Tube Size	Aux. Ging. Tube Size	Headgear Tube Size	Torque		Tip		Distal Offset		Left/Right	Weldable*
			Main	Aux.	Main	Aux.	Main	Aux.		
.022" x .028" x 4.0mm	.036" x 6.0mm	.045" x 4.5mm	-10°	—	0°	0°	12°	12°	L	391-094W
									R	391-093W

*Weldable tubes are only available individually welded to TPO bands.

Tip-Edge® System Molar Tubes

Slant Back with Easy-Out Feature

Combination


6/6 UPPER FIRST MOLAR LOW PROFILE

Rect. Tube Size	Round-Ging. Tube Size	Torque		Tip		Distal Offset		Left/Right	Weldable*	Bondable Big Foot Pad .325"
		Rect.	Rnd.	Rect.	Rnd.	Rect.	Rnd.			
.022" x .028"	.036"	-11°	—	0°	0°	6°	6°	L	391-872	391-907
x 4.7mm	x 6.0mm							R	391-873	391-908

Package of 10

*Weldable Tip-Edge tubes can be ordered individually on TPO bands; add W to part number.

Upper Right



6/6 LOWER FIRST MOLAR LOW PROFILE

Rect. Tube Size	Round-Ging. Tube Size	Torque		Tip		Distal Offset		Left/Right	Weldable*	Bondable Big Foot Pad .325"
		Rect.	Rnd.	Rect.	Rnd.	Rect.	Rnd.			
.022" x .028"	.036"	-27°	—	0°	0°	6°	6°	L	391-878	391-909
x 4.7mm	x 6.0mm							R	391-879	391-910

Package of 10

*Weldable Tip-Edge tubes can be ordered individually on TPO bands; add W to part number.

Lower Right



Combination Convertible

6/6 UPPER FIRST MOLAR

Rect. Tube Size	Round Tube Size	Torque		Tip		Distal Offset		Left/Right	Weldable*	Bondable	
		Rect.	Rnd.	Rect.	Rnd.	Rect.	Rnd.			Standard Pad	Big Foot Pad .325"
.022" x .028"	.036"	-10°	—	0°	0°	12°	12°	L	391-372	391-407	391-418
x 4.0mm	x 6.0mm							R	391-373	391-408	391-419

Package of 10

*Weldable Tip-Edge tubes can be ordered individually on TPO bands; add W to part number.

Upper Right



6/6 LOWER FIRST MOLAR

Rect. Tube Size	Round Tube Size	Torque		Tip		Distal Offset		Left/Right	Weldable*	Bondable	
		Rect.	Rnd.	Rect.	Rnd.	Rect.	Rnd.			Standard Pad	Big Foot Pad .325"
.022" x .028"	.036"	-27°	—	0°	0°	8°	8°	L	391-378	391-409	391-420
x 4.0mm	x 6.0mm							R	391-379	391-410	391-421

Package of 10


*Weldable Tip-Edge tubes can be ordered individually on TPO bands; add W to part number.

Lower Right



The Easy-Out feature straightens the archwire as it is pulled through the tube, lessening the force needed to remove the wire and improving patient comfort. Easy-Out also permits distal crown tipping on the mounted tooth.

Easy-Out Feature



Lingual Attachments

Bondable

Bondable Lingual Button

- For engaging elastomeric thread.
- Available on 2.16 x 3.8mm (.085" x .150") rectangular base for bicuspid and molars.
- Available on 3.6mm (.140") diameter base for cuspids.

224-112 Round Base (anterior)

224-113 Rectangular Base (posterior)
Package of 10



Bondable Plastic Button

- 3.1mm (.120") diameter head.

224-800 Plastic Lingual Buttons
Package of 10



Bondable Lingual Hook

224-116L Upper Left/Lower Right Bondable Hooks

224-116R Upper Right/Lower Left Bondable Hooks
Package of 10




Bondable Lingual Cleat

- Available with 3.6mm (.140") diameter base for cuspids.
- Available with 2.16 x 3.8mm (.085" x .150") rectangular base for bicuspid and molars.

224-114 Round Base (cuspids)

224-115 Rectangular Base (bicuspid and molars)
Package of 10



Lingual Attachments Weldable

One-Piece Lingual Button

- Head size 2.18mm (.086").

224-010 Flat Base
224-020 Curved Base
 Package of 10



Molar Cleat

- Stamped from 0.30mm (.012") stainless steel.
- Total length 7mm (.275").

224-004 Molar Cleats
 Package of 100



Cleat-Lug

- Formed from 0.30mm (.012") stainless steel.
- Length is 5mm (.200").

224-017 Posterior Cleat-Lug
 Package of 100



Ball-End Lingual Hook

226-013 Long Hook, Short Offset (A)
226-011 Short Hook, Short Offset (B)
226-033 Long Hook, Long Offset (C)
226-031 Short Hook, Long Offset (D)
 Package of 10



Lingual Cleat

224-007 Short, 4.4mm (.173")
224-008 Long, 5.2mm (.206")
 Package of 100



Lingual Molar Sheath

- 0.91mm (.036") I.D. for use with palatal arch bars.

328-054 Sheath with Latch Indent
328-254L Left Sheath with Distal Ball Hook
328-254R Right Sheath with Distal Ball Hook
 Package of 10



GripTite® First Molar Bands

- Mechanical interlocking indentations that are three times deeper than laser marked bands.
- Anatomically correct according to tooth shape and size.
- Available in 38 sizes,
- All bands are permanently identified with an advanced marking system that clearly indicates Palmer Notation, band size and tooth identification on the mesial surface.




Molar Band Kits

Introductory Kit	100 bands (per arch)
Standard Kit	200 bands (per arch)
Inventory Kit	400 bands (per arch)

Sizes 0 and 00 are not included in molar kits unless specifically requested.
Band boxes are furnished at no charge with all kit orders.

Upper and Lower Kit Distribution



GRIP TITE MOLAR BANDS – 38 SIZES

Upper/Lower Left/Right			.010 Firm Temper	.007 Standard Temper
First Molar	Upper	R	964-725	966-725
		L	964-726	966-726
	Lower	R	965-825	967-825
		L	965-826	967-826

Band Boxes

- Hinge top box with a divided plastic insert with each band kit order.
- Convenient band storage.

Upper Band Boxes—No Bands		Order No.
Bicuspid	Universal	100-955

Lower Band Boxes—No Bands		Order No.
1st Bicuspid	Universal	100-968
2nd Bicuspid	Universal	100-969

Pre-Fit® Second Molar Bands

- Anatomical right and left shapes.
- Available in 38 sizes (00 to 36).
- .005 mm thick.
- incorporate slight buccal indents and gingival interproximal festooning to promote positive seating at the proper levels.
- Advanced marking system clearly indicate the Palmer Notation, band size and tooth identification on the mesial.



Molar Band Kits

Introductory Kit	100 bands (per arch)
Standard Kit	200 bands (per arch)
Inventory Kit	400 bands (per arch)


Sizes 0 and 00 are not included in molar kits unless specifically requested.
Band boxes are furnished at no charge with all kit orders.

STANDARD TEMPER PLUS PRE-FIT MOLAR BANDS

	Upper/Lower	Left/Right	Order No.
2nd Molar	Upper	R	946-725
		L	946-726
	Lower	R	947-825
		L	947-826

Add standard band size numbers (00-36) after order numbers.

Upper and Lower Kit Distribution



Pre-Fit® Bicuspid Bands

- Markings indicate the Palmer Notation.
- Band size and tooth identification on the mesial.

Bicuspid Band Kits


Introductory Kit	150 bands (per arch)
Inventory Kit	300 bands (per arch)

UNIVERSAL BICUSPID BANDS


No Right and Left Shapes	Order No.
Upper First or Second	940-105
Lower First	942-107
Lower Second	942-105

There is no welding charge for custom-welding of TPO attachments.
A minimum order of 10 bands is requested, assortable by size and type.

Upper Kit Distribution



Lower Kit Distribution



Prewelded Band Assemblies

Many of our brackets, molar tubes and auxiliaries are available pre-welded to TPO's bands (we do not offer welding services for non-TPO parts).

Standard Welding Positions are used for fabricating all assemblies. GripTite Bands are used for first molars. Pre-Fit bands are used for second molars and bicuspid.


Standard Welding Positions

MOLARS

(Pre-adjusted, Standard Edgewise and Tip-Edge® prescriptions)


Upper Single Tube

Mesiodistal: Mesial end of tube is centered on mesial buccal cusp.
Long axis of tube is parallel to a line bisecting contact points.
Occlusogingival: Centered



Upper Combination Tube

Mesiodistal: Mesial end of tube is centered on mesial buccal cusp
Long axis of tube is parallel to a line bisecting contact points.
Occlusogingival: Rectangular tube centered. Round tube occlusal.



Lower Single Tube

Mesiodistal: Same as upper
Occlusogingival: Centered



BICUSPIDS

Upper

Mesiodistal: Centered
Occlusogingival: Offset occlusally



Preadjusted and
Standard Edgewise



Tip-Edge and
Tip-Edge PLUS

Lower

Mesiodistal: Centered
Occlusogingival: Offset gingivally



Note: The conical shape of Pre-Fit bicuspid bands gives them an "upside-down" appearance.

Prewelded Band Assemblies

Featuring GripTite Molar Bands

6 | 6 **GRIPTITE BAND WITH UPPER FIRST MOLAR MINI SIZE OCCLUSAL HEADGEAR TUBE**

Main Tube Size	Headgear Tube Size	Torque	Tip	Distal Offset		Left/ Right	Band Assemblies	
				Main	HG		Standard	Firm
.018" x .025" x 3.5mm	.045" x 3.5mm	-10°	0°	7°	0°	L		964-502
						R		964-503
.022" x .028" x 3.5mm	.045" x 63.5mm	-10°	0°	14°	0°	L	966-508*	964-508
						R	966-509*	964-509

Add band size number (00-36) after item number.

Tubes on page 43.

6 | 6 **GRIPTITE BAND WITH LOWER SINGLE FIRST MOLAR TUBE WITH DISTAL EXTENSION**

Archwire Tube Size	Tube Length	Torque	Tip	Distal Offset	Left/ Right	Band Assemblies	
						Standard	
.022" x .028"	4.2	-25°	0°	4°	L	967-628*	
					R	967-629*	

Add band size number (00-36) after item number.

Tubes on page 40.

6 | 6 **GRIPTITE BAND WITH LOWER SINGLE CONVERTIBLE FIRST MOLAR TUBE**

Archwire Tube Size	Tube Length	Torque	Tip	Distal Offset	Left/ Right	Band Assemblies	
						Firm	
.018" x .025"	4.2	-25°	0°	4°	L	965-536	
					R	965-537	
.022" x .028"	4.2	-25°	0°	4°	L	965-540	
					R	965-541	

Add band size number (00-36) after item number.

Tubes on page 41.

*Pre-welded band assemblies include lingual cleat.