



TP Orthodontics, Inc.

tportho.com

cosmeticbraces.com

ADHESIVES

ADHESIVES

Bonding Adhesive

- eXact® Orthodontic Adhesive

Banding Cement

- CrossLink® Light Cure Caplets
- CrossLink® Light Cure Syringes

Bonding Preparation

- TPO® Blue Etchant Gel
- Enamel Conditioner Solution
- Enamel Conditioner Gel
- Porcelain Primer

Adhesive Accessories

- Adhesive Mixing and Orientation Tray
- Tongue-Away® Shields
 - Clear Plastic
 - Colored Plastic
 - Disposable
- Lip Retractors
- FLEX-Retractor

eXact® Clear Orthodontic Adhesive

- Medium-high viscosity light cure adhesive.
- Use for metal, plastic or ceramic attachments.
- Cures completely clear - making it an exceptional system for aesthetic brackets.
- Prevents bracket drift, produces little flash and is extremely easy to clean up.



eXact Clear Light Cure Kit

151-270

Kit Contains:

- Light Cure Syringes, 4g (.141 oz), 2 ea.
- TPO Blue Etchant Gel in Syringe, 5 cc
- Replacement Tips (for Blue Etchant Gel)
- Light Cure Sealant, 6 cc
- Disposable Brushes, 50 ea.
- Foam Pellets, 100 ea.
- Spatula, 1 ea.
- Work Pad, 50 sheets

Replacement Materials

- 151-271** Paste in Syringe, 4g (.141 oz), 2 ea.
- 150-713** TPO Blue Etchant Gel in Syringe, 5 cc
- 151-272** Light Cure Sealant, 6 cc
- 150-704** Disposable Brushes, 50 ea.
- 150-705** Foam Pellets, 150 ea.
- 150-904** White Spatulas, 50 ea.
- 150-903** Work Pad, 50 sheets
- 100-620** Adhesive Mixing and Orientation Trays, 25 ea.

Single Syringe Package

- 151-273** Paste in Syringe, 2g (.0705 oz)

CrossLink® Light Cure Band Cement

Single Patient Caplets

- No-mix, no-mess, single-use dispensing caplets.
- Cement is blue to help clinicians with clean up.
- Disposable squeeze caplet has an angled, tapered stem for excellent control when dispensing paste.
- Cut the caplet tip and dispense cement into bands.
- Fluoride ion releasing paste does not require sealant, and there is no rush to beat setting time.
- Complete light cure is accomplished in 30 seconds.
- CrossLink provides excellent band retention.



Single Band Caplet

151-284 Small Patient Caplets, 30 caplets, .4g (.014 oz) each
Enough for one molar band

Preloaded Syringes

- Composite no-mix, no-mess cement, with fluoride ion release, preloaded in syringes.
- Cement is blue to help clinicians with clean up.
- No sealant is needed, and there is no rush to beat a setting time.
- Complete light cure is accomplished in 30 seconds.

151-285 CrossLink Kit - Blue, 24g (.847 oz) paste, 6 ea.



Replacement Materials

151-286 Blue Paste in Syringe, 4g (.141 oz)

150-904 White Spatulas, 50

Bonding Preparation

TPO® Blue Etchant Gel

- Can be used with any bonding system.
- Pre-loaded in a syringe.
- Curved tip simplifies application on high cuspids, bicuspids and molars.
- Disposable tips are easy to use, save time and eliminate the possibility of cross-contamination or infection.
- Etches enamel in 30 seconds.
- Phosphoric acid content 37%.



150-713 TPO Blue Etchant Gel in Syringe with Tips

Package of 1 syringe with 10 disposable tips

150-714 Disposable Dispensing Tips

Package of 25

Enamel Conditioner Solution

Medium Viscosity

- Can be used with any bonding system.
- Contains a plaque disclosing agent (green color) on the teeth after rinsing if plaque is still present.
- After cleaning plaque, patient should be re-pumiced and re-conditioned.

150-020 Enamel Conditioner Solution, 9g (.317 oz)

Phosphoric acid content 37%



Enamel Conditioner Gel

Thick Viscosity

- Can be used with any bonding system.
- Thicker than enamel conditioner solution.
- Will not overflow onto gingival tissue when conditioning enamel.
- Blue color is readily seen during application.
- Easily removed by evacuation and/or rinsing with water.

150-015 Enamel Conditioner Gel, 20g (.705 oz)

Phosphoric acid content 37%



Porcelain Primer

- Prepares the surface of glazed dental porcelain for bonding of orthodontic attachments.
- The silane coupling agent chemically unites to the bonding adhesive and the silicon in porcelain.

150-810 Porcelain Primer, 16 cc



Adhesive Accessories

Adhesive Mixing and Orientation Tray

- Organizes bondable brackets and tubes for an individual patient.
- Adequate area available to mix bonding adhesive.
- Small integrated cups hold conditioner or paste components.

100-620 Adhesive Mixing and Orientation Tray
Package of 25



Tongue-Away® Shields

- Restricts the tongue to prevent moisture contamination and alleviates patient fatigue.
- Vertical shield is offset to provide adequate tongue space.
- Surfaces are angulated and serrated for retention.
- Holes facilitate placement and removal with a bird-beak plier, permit attachment of a safety string or allow insertion of a saliva ejector.

Clear Plastic Tongue-Away Shields

100-398 Tongue-Away Kit
4 shields, one of each size

100-395 Small

100-396 Medium

100-397 Large

100-496 Super
Package of 3



Phthalate-Free

Colored Plastic Tongue-Away Shields

100-390 Blue

100-391 Green

100-392 Red

100-393 Blue, Green, Red
Package of 3, medium size only



Phthalate-Free

Disposable Tongue-Away Shields

100-494 Polystyrene
Package of 10, medium size only



Adhesive Accessories

Lip Retractors

- Provide view of the full arch for intraoral photography or when etching, direct bonding or banding.
- Holes assist patients to firmly hold retractors.
- Made from cold sterilizable plastic.

100-499 Adult Size Retractors Only

100-498 Child Size Retractors Only
One Pair



FLEX-Retractor

- One-piece design.
- Clear, easy-to-maintain shatterproof plastic.
- Fit most patients.
- Easy to insert and manipulate.
- Provide a view of the full arch for direct bonding procedures or intraoral photography.

100-517 FLEX-Retractor
Package of 3



