








Wire

The products covered by this instruction include:

- Nickel Titanium Archwire
- Stainless Steel Archwire
- Titanium Archwire
- Braided Wire,
- Tip-Edge System Archwire
- Ligature Wire
- Specialty Archwires

Symbols Used on Labeling

	Manufacturer.
	Date of manufacture.
	Batch code.
	Caution! See warnings and precautions.
	Do not re-use.
	Use by date.
	Do not use if package is damaged.

Indications for Use

Intended for use in orthodontic treatment.

Warnings

- This product contains nickel and/or chromium which has been known to cause allergic reactions.
- This is a single-use device. Any reuse of this product will run the risk of cross-contamination to patients.

Precautions

- This product must be used by qualified orthodontists or dentists trained in orthodontics.

Contraindications

- Installing the device to patients with compromised oral health, such as poor oral hygiene, periodontal disease, or preexisting dental conditions (i.e., decalcification, carious lesions, large restorations) can increase the risk of tooth damage.
- TP Orthodontics sells its products to trained dental professionals and orthodontists, who have the primary responsibility to identify any possible contraindications that may preclude the use of this product.

Instructions for Use

This product must be used by qualified orthodontists or dentists trained in orthodontics, who are knowledgeable in its use.

Storage Condition / Handling

Product should be kept in a clean location.

Sterilization of these single use devices is not required. The clinician (qualified orthodontists or dentists trained in orthodontics) is responsible for assessing the specific needs for their patients / practice.

Product Performance Structure and Composition

The wires are composed of either stainless steel, nickel titanium or titanium alloys.

Aesthetic wires additionally contain a teflon coating.

NOTICE:

TP Orthodontics, Inc. warrants that its products are free from defects in materials and workmanship. TPO's liability and the customer's sole remedy in the event of any claimed defect shall be limited to replacement of the product, or at TPO's discretion, refund of the purchase price. This warranty is in lieu of any other warranties, express or implied, including any implied warranty of merchantability, fitness for a particular purpose, or otherwise, and of any other obligations or liability on TPO's part. Other than the warranties listed above, under no circumstances will TPO be liable for any loss, damages arising in connection with use, or inability to use, TPO's products.

