

Nu-Edge® SL

Brackets de cromo-cobalto



Directrizes de sequência de arcos

1

ALINHAMENTO

.014"

TPO® Original Reflex® Nickel Titanium Wire
TPO Reflex Heat-Activated Nickel Titanium Wire

.016"

TPO Original Reflex Nickel Titanium Wire
TPO Reflex Heat-Activated Nickel Titanium Wire

2

NIVELAMENTO

.016" x .022"

TPO Original Reflex Nickel Titanium Wire
TPO Reflex Heat-Activated Nickel Titanium Wire
TPO TiMolium® Titanium Wire

.018" x .025"

TPO Original Reflex Nickel Titanium Wire
TPO Reflex Heat-Activated Nickel Titanium Wire

.017" x .025"

TPO TiMolium Titanium Wire

3

FASE DE TRABALHO

.016" x .022"

TPO Shiny Bright Stainless Steel Wire

.018" x .025"

TPO Shiny Bright Stainless Steel Wire

4

DETALHES DA FASE DE ACABAMENTO

.016" x .022"

TPO Shiny Bright Stainless Steel Wire

.018" x .025"

TPO Shiny Bright Stainless Steel Wire

TOTAL ENGAJAMENTO DO SLOT

.017" x .025"

TPO Shiny Bright Stainless Steel Wire

.0215" x .025"

TPO Shiny Bright Stainless Steel Wire

Mais informações em **tportho.com**
ou ligue **+351 21 124 6897**

