

















# INSTRUCTIONS

## Labeling Symbol Glossary for TP Orthodontics

TP Orthodontics utilizes symbols in compliance with ISO 15223-1 Medical Devices – Symbols to be supplied by the manufacturer. Additional symbols added to comply with regulatory requirements are listed below. These symbols are found in the product packaging, labels and/or instructions for use.

	Manufacturer.		Date of manufacture.
	Catalog number		Batch code.
	CE Marking – indicates device complies to EU regulatory requirements	<b>Qty.</b>	Quantity
<b>0297</b>	Notified Body Number – Accompanies CE mark of class IIa devices		Authorized Representative in the European community.
	Medical device.	<b>Rx ONLY</b>	Prescription device.
	Caution! See warnings and precautions.		Consult instructions for use.
	Do not re-use.		Use by date.
	Keep away from sunlight.		Store between 2°C/35°F (lower limit) and 27°C/80°F (higher limit).
	Do not use if package is damaged.		Contains hazardous substance – indicates device contains cobalt >0.1% of total weight
	Contains or presence of natural rubber latex		